



Recovery activated



2025

Instruction Videos



PlasmaFlow®



PlasmaWave®

Your patients simply scan the QR code with a smartphone to view a video about the ManaMed product you prescribed. The videos show how to apply and use the product, and also include tips and troubleshooting.



Kangaroo Sling



Kahuna Brace



ManaFlow® 51



ManaFlow® 52



ManaFlexx 2®



ManaEZ® ROM Ice



ManaEZ® Hip Ice



Mini Kahuna Brace



Bionic Sling



GOAT



Tailback Line



Tailback 64



ManaEZ® Shoulder Ice



ManaEZ® ROM Elite/
ManaEZ® ROM Elite Short



ManaEZ® Tendon Strap



ManaEZ® Humerus Uni



ManaEZ® Boot Air



ManaEZ® Collar



ManaEZ® Soft Collar



ManaEZ® Spring &
Standard Crutch



ManaEZ® Hinged Knee



ManaEZ® Thumb Spica



Royal Boot



ManaEZ® Ankle Stirrup



ManaEZ® Ankle 02



PDAC Letters

Disclaimer: The HCPCS codes listed in this catalog are provided as suggested codes based on publicly available information. We do not make claims, promises, or guarantees as to the availability of reimbursement for any product. It is within the sole discretion of the customer to determine the appropriate billing code for a product, as well as whether the use of a product complies with medical necessity and other documentation requirements of the payor. This information is not intended to be, nor should it be considered, billing or legal advice.

The Company

Welcome to the new chapter of ManaMed, where pioneering orthopedic solutions meet the spirit of "Recovery Activated," driving us to the forefront of patient care and mobility. Our comprehensive product portfolio, spanning vascular therapies, assisted mobility, and DVT prevention, is designed to streamline the rehabilitation process for those eager to reclaim their active lifestyles. Our unwavering commitment to innovation keeps us constantly exploring, learning, and pushing the boundaries of what's possible in orthopedic care.

Central to our achievements are the significant partnerships we've nurtured within the industry, enhancing the level of care and offering products that expedite the return to normalcy for patients. Highlighting this is our esteemed collaboration with the UFC Performance Institute, ensuring that UFC athletes have access to the highest quality orthopedic braces, soft goods, and essential stabilizing equipment like crutches, knee braces, and walking boots. Similarly, our dedication to our distributor network—a robust community committed to excellence—reflects our ongoing investment in manufacturing, quality control, and customer service excellence.

Our rebranding transcends the visual aspects to embody our deepened pledge to not only meet but exceed the expectations of our clients and their patients. As ManaMed continues to expand, we remain your steadfast ally in orthopedic recovery, always ready to face new challenges with creative solutions. This journey is not just about us; it's about fostering enduring connections with our customers, who are the heart and soul of our company. Our success is deeply rooted in our ability to truly understand and satisfy our customers' needs. When you partner with ManaMed, you become an integral part of a mission dedicated to promoting healthier, happier lives.



At ManaMed, we believe that our customers are the heart and soul of our company. Our mission is to design, develop, and distribute orthopedic solutions that keep active, high-performing individuals in the game and back on track following injury or surgery. Our extensive catalog includes vascular and electrotherapies, assisted mobility, electrical stimulation, and DVT prevention products. By partnering with orthopedists, physical therapists, surgeons, and DME companies across the US, we ensure that patients of all ages and backgrounds remain healthy and comfortable on their path to recovery. When you choose us, you're not just a customer – you're part of our mission to help people live healthier, happier lives.

Trevor Theriot

President and CEO

Online Ordering Portal

Our one stop shop for ordering, paying invoices, and tracking orders.

View your order history and create ordering templates for easy reordering. Quickly pay invoices via ACH, credit card, or apply your credit memos.

Contact us at manamed.com/become-a-customer and we will get you set up.

Custom Logo Program

Add your logo to the products you order on our portal to create better brand visibility. We keep an inventory of your logos and work with you to keep them in stock.

We use stickers for hard goods and heat transfers for soft goods in three different sizes to match our various products.

Available for most products, contact us at logoprogram@manamed.com for more information.



Table of Contents

Vascular Therapy	PlasmaFlow	6	Foot and Ankle	MAC Boot	19
	PlasmaFlight	6		Royal Boot Air/Royal Boot Air Short	19
	PlasmaWave	6		ManaEZ Boot Air/ManaEZ Boot Air Short	19
	ManaEZ Compression Sock	6		ManaEZ Boot Pediatric	20
				ManaGO 71	20
				ManaEZ Ankle 71	20
				ManaEZ Ankle 06	20
				ManaEZ Ankle 02	21
				ManaEZ Ankle Stirrup	21
				ManaEZ Post-Op SQ Toe Shoe	21
				ManaEZ Night Splint	21
Lymphedema Therapy	ManaFlow 51	7			
	ManaFlow 51 Full	7			
	ManaFlow 52	7			
	ManaFlow 52 Full	7			
Surgical Bundles and Supplies	Total Joint Bundle	9	Spine	Tailback 2-64	23
	Spine Bundle	9		Tailback 2-56	23
	Ankle Bundle	9		Tailback 2-50	23
	Hip Bundle	9		Tailback 2-48	24
	Shoulder Bundle	10		Tailback 2-42	24
	Cervical Bundle	10		Tailback 27	24
	ManaMomma Maternity Bundle	10		Halfback 50	24
	Reddy Joint Bath	10		ManaMomma Maternity Back Brace	25
				ManaEZ Collar 80	25
				ManaEZ Collar 74	25
				ManaEZ Soft Collar	25
Shoulder	Kahuna Brace	12	Pain and Recovery Solutions	ManaSport+	27
	Mini Kahuna Brace	12		ManaFlexx2	27
	ManaEZ Sling	12			
	ManaEZ Humerus	13	Circulating Cold Therapy Devices	ManaCold Platinum 2.0 Universal	28
	Mana Water Sling	13			
Elbow	Cabo Elbow	13	Braces with Ice and Compression	ManaEZ Hip Ice	29
				ManaEZ ROM Ice	29
Hand and Wrist	ManoFX	14		ManaEZ Hinged Ice	29
	ManaEZ Thumb Spica	14		ManaEZ Ankle Ice	29
	ManaEZ Wrist	14		ManaEZ Shoulder Ice	30
				ManaEZ Ice Wrist	30
Knee	ManaEZ ROM Elite/ManaEZ ROM Elite Short	16	Ice Wraps	ManaEZ Ice Wrap Hip	30
	ManaEZ Hinged Knee	16		ManaEZ Ice Wrap Knee	30
	ManaEZ Wrap 33	16	Mobility Aids	ManaEZ Crutch	31
	ManaEZ Knee Immobilizer	16		ManaEZ Crutch Spring	31
	ManaEZ Tendon Strap Universal	17	Maternity	ManaMomma Maternity Back Brace	31
	EZOA	17		Three-In-One Belly Band	31
	GOAT	17			
	Cabo Knee	17			

Vascular Therapy





PlasmaFlow®

PlasmaFlow®

The PlasmaFlow is intended to be an easy to use sequential compression system prescribed by a physician. For use in the home or clinical setting to help prevent the onset of DVT in patients, by stimulating blood flow in the extremities by simulating muscle contractions. This device is used to aid in preventing DVT by enhancing blood circulation, diminishing post-operative pain and swelling, and enhancing wound healing.

Features:

- Long-lasting rechargeable battery. Up to 10 hours of use on one charge
- 2 Modes: Continuous inflation and Step-Up inflation technology which allows the unit to increase pressure in slower increments to 55mmHG
- 2 LCD screens to monitor usage and pressure
- Portable, lightweight, and tubeless
- Suggested HCPCS: E0676

Product Name	Part Number
PlasmaFlow DVT	PF0001
PlasmaFlow Extension Sleeves	PFEXT
PlasmaFlow Charger	PFCHG
PlasmaFlow Hospital Charger	PFHCHG



PlasmaFlight®

PlasmaFlight®

The PlasmaFlight® has been approved by the FDA as an over-the-counter portable inflatable massage system that simulates the kneading and stroking of tissue with the hands by use of an inflatable cuff. This device is used for increasing blood circulation in the treated areas, as well as relief of minor muscle aches and pains.

Features:

- Pump will inflate each leg cuff to pre-set pressure of 55mmHG and releases once that pressure is reached
- Visual display shows hourly usage, mode, and pressure level
- Electronically controlled pumps compress the calf or calves to aid in blood flow in the lower extremity
- Internal rechargeable battery and cordless design allow for ease of placement and patient mobility
- Suggested HCPCS: E0676

Product Name	Part Number
PlasmaFlight	PFLIGHT01



PlasmaWave®

PlasmaWave®

The PlasmaWave is a game-changing clinical grade pneumatic SCD device used in the home or clinical setting. It has a market-leading 25+ hour battery.

Features:

- Over 8 feet of reusable tubing
- 25+ hour battery
- Sleeve options include both reusable or disposable
- Convenient bed rack also available
- 4 Unique pressure modes which allows the unit to increase pressure in slower increments
- 1 LCD screen to monitor usage, pressure, and battery life
- Lightweight
- Suggested HCPCS: E0676

Product Name	Part Number
PlasmaWave Portable DVT	PWAVE01
PlasmaWave Bed Rack	PWRACK
PlasmaWave Bed Holder	PWBH01
PlasmaWave Charger	PWCHG
PlasmaWave Tube Set	PW15
PlasmaWave Standard Calf Sleeves	M-40
PlasmaWave XL Calf Sleeves	M-42
PlasmaWave Standard Foot Garment	M-21
PlasmaWave Standard Thigh Garment	M-50
Go Green Sleeve Barrier (10 pairs per box)	M-32
M-31 Reusable Universal Calf Garment	M-31

M-31 Reusable Sleeve
• Suggested HCPCS: E0673



Go Green
Compression Sleeve Barrier



ManaEZ® Compression Sock

ManaEZ Compression Sock

The ManaEZ Compression sock is a quick fitting, sized sock that uses a gradient descent compression design. The sock is made with breathable nylon/spandex material that is quick-drying and has enhanced air circulation and odor fighting properties.

- Latex free
- Unique 3D stereoscopic pattern design that enhances shock absorption
- Machine washable with cold water
- Suggested HCPCS: A6530

Product Name	Part Number
ManaEZ Compression Sock Medium	EZSOCK01M
ManaEZ Compression Sock Large	EZSOCK01L
ManaEZ Compression Sock Extra Large	EZSOCK01XL



ManaFlow® 51

ManaFlow® 51 (Uncalibrated)

The ManaFlow 51 provides a preset 50 mmHG in each of the four chambers in our tubeless pneumatic compression device used to treat lymphedema or aid in the prevention of DVT.

Features:

- Rechargeable battery takes 3 hours to charge, 8-hour use on a single charge, can be used while plugged in
- Compliance meter shows use, pressure, and chamber
- Gradient compression
- Easy-to-clean material
- Suggested HCPCS: E0651 and E0669

Product Name	Part Number
ManaFlow 51 (Uncalibrated)	MFLOW51
ManaFlow 51/52 Extenders	MFLOWEXT



ManaFlow® 51 Full

ManaFlow® 51 Full (Uncalibrated)

The ManaFlow 51 Full provides a preset 50 mmHG in each of the four-chambers in our tubeless pneumatic compression device used to treat lymphedema or aid in the prevention of DVT. Our full version extends from the top of the thigh to the foot.

Features:

- Rechargeable battery takes 3 hours to charge, 8-hour use on a single charge, can be used while plugged in
- Compliance meter shows use, pressure, and chamber
- Gradient compression
- Easy-to-clean material
- Suggested HCPCS: E0651 and E0667

Product Name	Part Number
ManaFlow 51 Full (Uncalibrated)	MFLOW51LEG



ManaFlow® 52

ManaFlow® 52 (Calibrated)

The ManaFlow 52 allows the clinician to calibrate the pressure of our four-chamber tubeless pneumatic compression device, used to treat lymphedema or aid in the prevention of DVT.

Features:

- Rechargeable battery takes 3 hours to charge, 8-hour use on a single charge, can be used while plugged in
- Compliance meter shows use, pressure, and chamber
- Gradient compression
- Easy-to-clean material
- Suggested HCPCS: E0652 and E0669

Product Name	Part Number
ManaFlow 52 (Calibrated)	MFLOW52
ManaFlow 51/52 Extenders	MFLOWEXT



ManaFlow® 52 Full

ManaFlow® 52 Full (Calibrated)


The ManaFlow 52 Full allows the clinician to calibrate the pressure of our four-chamber tubeless pneumatic compression device, used to treat lymphedema or aid in the prevention of DVT. Our full version extends from the top of the thigh to the foot.

Features:

- Rechargeable battery takes 3 hours to charge, 8-hour use on a single charge, can be used while plugged in
- Compliance meter shows use, pressure, and chamber
- Gradient compression
- Easy-to-clean material
- Suggested HCPCS: E0652 and E0667

Product Name	Part Number
ManaFlow 52 Full (Calibrated)	MFLOW52LEG

Surgical Bundles and Supplies



Surgical Bundles



Total Joint Bundle

Total Joint Bundle

This Total Joint Bundle is an industry first; this game-changing Bundle has everything your patient needs to recover from a surgery or injury.

Bundle contains:

- ManaEZ ROM Ice our cold compression device for pain, swelling, and range-of-motion control
- PlasmaFlow our Portable DVT prevention device
- ManaFlexx 2 our NMES/TENS device to prevent muscular atrophy

- Suggested HCPCS: L1833, E0676, E0745

Product Name	Part Number
Total Joint Bundle	TJK001



Spine Bundle

Spine Bundle

The Spine Bundle is a game changer for the DME or surgeon to provide everything your patient needs for an improved post-op outcome.

Bundle contains:

- Tailback 2-50 our LSO support
- PlasmaFlow our portable DVT prevention device
- ManaFlexx 2 our NMES/TENS device to help prevent muscular atrophy

- Suggested HCPCS: L0650, E0676, E0745

Product Name	Part Number
Spine Bundle	SK001



Ankle Bundle

Ankle Bundle

The Ankle Bundle is a game changer for the DME or surgeon to provide everything your patient needs for an improved post-op outcome.

Bundle contains:

- MAC Boot our pneumatic walker
- ManaEZ Ankle Ice our wrap-around ankle orthosis that has ice and compression for edema and pain relief
- PlasmaFlow our portable DVT prevention device
- ManaFlexx 2 our NMES/TENS device to help prevent muscular atrophy

- Suggested HCPCS: L4361, L1902, E0676, E0745

Product Name	Part Number
Ankle Bundle	AK001



Hip Bundle

Hip Bundle

The Hip Bundle is a game changer for the DME or surgeon to provide everything your patient needs for an improved post-op outcome.

Bundle contains:

- ManaEZ ROM Hip Ice our range-of-motion cold compression wrap-around orthosis for edema and pain control
- PlasmaFlow our portable DVT prevention device
- ManaFlexx 2 our NMES/TENS device to help prevent muscular atrophy

- Suggested HCPCS: L1686, E0676, E0745

Product Name	Part Number
Hip Bundle	HK001

Surgical Bundles *[continued]*



Shoulder Bundle

Shoulder Bundle

The Shoulder Bundle is a game changer for the DME or surgeon to provide everything your patient needs for an improved post-op outcome.

Bundle contains:

- Kahuna Brace our adjustable airplane shoulder support
- ManaEZ Shoulder Ice our wrap-around cold compression device for edema and pain control
- PlasmaFlow our portable DVT prevention device
- ManaFlexx 2 our NMES/TENS device to help prevent muscular atrophy
- Suggested HCPCS: L3960, L3670, E0676, E0745

Product Name	Part Number
Shoulder Bundle	SHK001



Cervical Bundle

Cervical Bundle

The Cervical Bundle is a game changer for the DME or surgeon to provide everything your patient needs for an improved post-op outcome.

Bundle contains:

- ManaEZ Collar 80 our adjustable multi-post collar with mandibular and occipital support
- PlasmaFlow our portable DVT prevention device
- ManaFlexx 2 our NMES/TENS device to help prevent muscular atrophy and provide pain relief
- Suggested HCPCS: L0180, E0676, E0745

Product Name	Part Number
Cervical Bundle	CK001



ManaMomma Maternity Bundle

ManaMomma Maternity Bundle

The ManaMomma Maternity Bundle is a game changer for the DME or provider to supply everything your patient needs for an improved outcome pre and postpartum.

Bundle contains:

- ManaMomma Maternity Back Brace – spinal orthosis for pendulous abdomen
- Three-In-One Belly Band – elastic abdominal and pelvic support system
- ManaFlexx 2 – NMES/TENS unit with a pair of lead wires and electrodes
- ManaEZ Compression Socks – gradient compression sock
- Suggested HCPCS: E0745, A6530, L0650, L0621

Product Name	Part Number
ManaMomma Maternity Bundle	MMB001

Surgical Supplies



Reddy Joint Bath

Reddy Joint Bath

CDC-approved method of wound lavage with our povidone-iodine solution for PJI and SSI

Features:

- Industry-leading wound lavage
- 3-minute soak time
- Certified sterile
- Bundle includes: 5 saline solution bottles, iodine packet, bowl, and prep tray

Product Name	Part Number
Reddy Joint Bath	MMJB001

Upper Extremity





Kahuna Brace

Kahuna® Brace

The Kahuna Brace is a universally-sized abduction sling and is applied to the left or right side.

Features:

- Multiple positioning options; the Kahuna malleable aluminum waistband with arm support allows for unlimited arm positioning including airplane splint
- Breathable material allows for optimal comfort
- Ideal for shoulder injuries and their positioning needs
- Universal sizing reduces your inventory
- Foam ball for enhanced circulation
- Step-down abduction pillow available for purchase

- Suggested HCPCS: OTS L3999 or CF L3960

Product Name	Part Number
Kahuna Brace	KB001
Replacement Kahuna Pillow	KB001-P
Replacement Kangaroo Pillow	KSS001-P



Mini Kahuna® Brace

Mini Kahuna® Brace

The Mini Kahuna Brace is a universally-sized abduction sling and is applied to the left or right side.

Features:

- The Kahuna Pillow included is a standard black abduction pillow
- Pillow conforms to the body and is secured with a waist strap
- Breathable material allows for optimal comfort
- Ideal for shoulder injuries and their positioning needs
- Universal sizing reduces your inventory
- Foam ball for enhanced circulation

- Suggested HCPCS: OTS L3999 or CF L3670

Product Name	Part Number
Mini Kahuna Brace	MKB001
Replacement Kahuna Pillow	KB001-P
Replacement Kangaroo Pillow	KSS001-P



ManaEZ® Sling

ManaEZ® Sling

The ManaEZ Sling is a shoulder immobilizer and includes a waist strap.

Features:

- Soft low-profile cuff
- Soft padded shoulder strap for added comfort
- Waist strap ensures proper arm position
- Is used for mild sprains, strains, stable humeral fractures, and elbow injuries
- Sizes from Small to XXL
- Also available as a Universal size

- Suggested HCPCS: L3670

Product Name	Part Number
ManaEZ Sling Small (12.4" x 6.3")	EZSL01S
ManaEZ Sling Medium (14" x 7")	EZSL01M
ManaEZ Sling Large (15.5" x 7.9")	EZSL01L
ManaEZ Sling Extra Large (17" x 8.6")	EZSL01XL
ManaEZ Sling Extra Extra Large (18.5" x 9.5")	EZSL01XXL
ManaEZ Sling Universal	EZSLUNI

Shoulder *(continued)*



ManaEZ® Humerus

ManaEZ® Humerus Universal

The ManaEZ Humerus Universal is used to stabilize humeral fractures while allowing protected motion of the arm. Features include a shoulder cap design, customizable lightweight plastic with breathable foam, and a universal left or right arm with adjustable straps that fit most arms.

Features:

- Shoulder cap design
- Lightweight and breathable
- Mold-able and trim-able
- Universal left or right arm
- Adjustable straps fit any size
- Closed-cell foam for easy cleaning

- Suggested HCPCS: L3980

Product Name	Part Number
ManaEZ Humerus Universal	EZH3980



ManaWater Sling

ManaWater Sling

The ManaWater Sling is a waterproof arm sling that is used for maintaining arm position while bathing. Features include one-handed adjustment, universal left/right fit, and an adjustable wrist strap that fits any size.

Features:

- One-handed adjustment
- Lightweight and waterproof
- Adjustable wrist strap fits any size
- Can be worn on both right or left arms

- Suggested HCPCS: L9900

Product Name	Part Number
ManaWater Sling	MWSLING01

Elbow



Cabo® Elbow

Cabo® Elbow

This universal-fit range-of-motion elbow brace is a great alternative to casting or repetitive splinting. It treats any injury in which a limited range of elbow motion is desired.

Features:

- ROM hinge features flexion setting from 0-120 degrees and extension settings from 0-90 degrees
- Telescoping length
- Drop lock with three preset locking positions (0, 15, and 30 degrees)
- Lightweight design for comfort and ease of use
- Optional hand grip
- Adjustable straps for universal fit

- Suggested HCPCS: OTS L3761 or CF L3760

Product Name	Part Number
Cabo Elbow Universal Left	CE001L
Cabo Elbow Universal Right	CE001R

Hand and Wrist



ManoFX

ManoFX

The ManoFX is a wrist, hand, and finger orthosis. The unique design allows fractures of the hand and finger to be positioned for optimum healing. This innovative approach allows customization of the brace to be configured as a ulnar gutter splint (boxer's fracture), radial gutter splint, or a metacarpal fracture brace (buddy finger splint).

Features:

- Malleable aluminum stays
- Split dorsal stay
- Soft breathable body
- Wraparound design for customized fit
- Suggested HCPCS: OTS L3809 or CF L3807

Product Name	Part Number
ManoFX Universal Size - Right	MANOFX-R
ManoFX Universal Size - Left	MANOFX-L



ManaEZ[®] Thumb Spica

ManaEZ[®] Thumb Spica

The ManaEZ Thumb Spica is a universal orthosis that includes a double speed-lacer closure.

Features:

- Low profile, high-quality, breathable material
- Three malleable thumb stays
- Palmer stay
- Lateral stay
- Radial stay
- Universally sized for left or right
- Suggested HCPCS: OTS L3809 or CF L3807

Product Name	Part Number
ManaEZ Thumb Spica Universal Left	EZTS3809L
ManaEZ Thumb Spica Universal Right	EZTS3809R



ManaEZ[®] Wrist

ManaEZ[®] Wrist

The ManaEZ Wrist is a high quality, universal wrist orthosis. Single speed-lacer closure makes it easier for the user to self apply the brace.

Features:

- Low profile, high-quality, breathable material
- Palmer stay
- Dorsal stay
- Universally sized for left or right
- Suggested HCPCS: L3908

Product Name	Part Number
ManaEZ Wrist Universal Left	EZW3908L
ManaEZ Wrist Universal Right	EZW3908R

Knee





ManaEZ® ROM Elite Short

ManaEZ® ROM Elite

ManaEZ® ROM Elite/ManaEZ ROM Elite Short

Features:

- ManaEZ ROM Elite has telescoping and malleable hinge bars extending from 16.5"–27" (42cm–68.5 cm)
- ManaEZ ROM Elite Short version offers a fixed length of 18"
- Hinge range of motion adjustment: Flexion (0° to 120°) and extension (0° to 90°)
- Drop lock offers the ability to control and lock the hinge motion to 0°, 15°, or 30°
- Lightweight, latex-free, and breathable material
- Easy application, universal size with quick-release buckles, and fits either left or right leg

- Suggested HCPCS:
OTS L1833 or CF L1832

Product Name	Part Number
ManaEZ ROM Elite (Size Universal)	EZRE001
ManaEZ ROM Elite Short (Size Universal)	EZRS1833



ManaEZ® Hinged Knee

ManaEZ® Hinged Knee

The ManaEZ Hinged Knee is a convenient wrap-on hinged sleeve. The open popliteal feature allows for a comfortable fit.

Features:

- Internal hinges enhance medial and lateral support
- Condyle pads for improved fit and added comfort
- Wrap feature for easy on/off application
- High-quality Neoprene
- Open Popliteal

- Suggested HCPCS:
OTS L2999 or CF L1820

Product Size	Thigh Circumference	Knee Circumference	Calf Circumference
S	15" – 18.5"	13" – 14"	12" – 14"
M	18.5" – 21"	14" – 15"	14" – 16"
L	21" – 23.5"	15" – 17"	16" – 18"
XL	23.5" – 26.5"	17" – 19"	18" – 20"
XXL	26.5" – 29.5"	19" – 21"	20" – 22"

Product Name	Part Number
ManaEZ Hinged Knee Small	EZHK1820S
ManaEZ Hinged Knee Medium	EZHK1820M
ManaEZ Hinged Knee Large	EZHK1820L
ManaEZ Hinged Knee Extra Large	EZHK1820XL
ManaEZ Hinged Knee Extra Extra Large	EZHK1820XXL



ManaEZ® Wrap 33

ManaEZ® Wrap 33

The ManaEZ Wrap 33 has unique anterior and posterior body closures, which allow the user to maintain proper hinge alignment while fitting.

Features:

- Universal sizing for left or right knee, reducing inventory
- Range of motion hinges for controlled motion
- Medial and lateral hinge provides excellent support
- Breathable neoprene body provides warmth and compression

- Suggested HCPCS:
OTS L1833 or CF L1832

Product Name	Part Number
ManaEZ Wrap 33 Universal	EZW33



Single Panel

3 Panel

ManaEZ® Knee Immobilizer

ManaEZ® Knee Immobilizer

The soft, breathable design accommodates most leg circumferences, both styles can be adjusted and trimmed to provide a customized fit. The combination of the posterior stay and two adjustable lateral stays maintains the optimum leg position. The popliteal pad ensures an intimate fit around the knee joint. Available in 18", 20", and 22" lengths in our single and three-panel design. Ideal for strains and sprains of the knee, stable fractures of the tibia, and patellar tendon repairs.

- Suggested HCPCS: L1830

Product Name	Part Number
ManaEZ Single Panel Knee Immobilizer 18"	EZIMM01-18
ManaEZ Single Panel Knee Immobilizer 20"	EZIMM01-20
ManaEZ Single Panel Knee Immobilizer 24"	EZIMM01-24
ManaEZ 3 Panel Knee Immobilizer 18"	EZ3PIMM01-18
ManaEZ 3 Panel Knee Immobilizer 20"	EZ3PIMM01-20
ManaEZ 3 Panel Knee Immobilizer 24"	EZ3PIMM01-24



ManaEZ® Tendon Strap

ManaEZ® Tendon Strap Universal

The ManaEZ Tendon Strap Universal is used to reduce tension on the patellar tendon.

- Low profile
- Neoprene strap
- Gel pad for compression
- Suggested HCPCS: L1825

Product Name	Part Number
ManaEZ Tendon Strap Universal	EZTS01



EZOA

EZOA

The EZOA is a lightweight, low-profile, single-upright OA orthosis, ideal for moderate to severe OA.

- Adjustable varus/valgus for customized pain relief
- Butterfly adjuster allows for easy fitting without tools
- Non-slip liner designed for all-day support
- Universally sized but left or right specific
- Malleable cuffs for optimum support
- Suggested HCPCS: OTS L1851 or CF L1843

Product Name	Part Number
EZOA Left Knee	OFF1L
EZOA Right Knee	OFF1R



GOAT

GOAT

The GOAT is a universally sized, low-profile, lightweight option for patients with mild to moderate knee OA.

- Butterfly varus/valgus adjustment, no tools needed
- One brace does medial and lateral offloading
- One size fits all
- Right and left specific
- Ideal for mild to moderate OA
- Breathable and low-profile neoprene
- Suggested HCPCS: OTS L1851 or CF L1843

Product Name	Part Number
GOAT Left	GOATL
GOAT Right	GOATR



CaboKnee®

Cabo® Knee

The Cabo Knee is a dual-upright knee brace, used for support and protection during high-demand activities.

- Lightweight aluminum alloy
- Quick-fitting buckles
- Comfortable non-slip liner
- Suggested HCPCS: OTS L1852 or CF L1845

Product Size	Thigh Circumference	Knee Circumference	Calf Circumference
S	15" - 18.5"	13" - 14"	12" - 14"
M	18.5" - 21"	14" - 15"	14" - 16"
L	21" - 23.5"	15" - 17"	16" - 18"
XL	23.5" - 26.5"	17" - 19"	18" - 20"
XXL	26.5" - 29.5"	19" - 21"	20" - 22"

Product Name	Part Number
Cabo Knee Right Knee Brace Small	TRBL45S1
Cabo Knee Left Knee Brace Small	TRBL45S2
Cabo Knee Right Knee Brace Medium	TRBL45M1
Cabo Knee Left Knee Brace Medium	TRBL45M2
Cabo Knee Right Knee Brace Large	TRBL45L1
Cabo Knee Left Knee Brace Large	TRBL45L2
Cabo Knee Right Knee Brace Extra Large	TRBL45XL1
Cabo Knee Left Knee Brace Extra Large	TRBL45XL2
Cabo Knee Right Knee Brace XXL Large	TRBL45XXL1
Cabo Knee Left Knee Brace XXL Large	TRBL45XXL2

Foot and Ankle





MAC – ManaMed Air Customizable Walking Boot

Built-in bladder system and release valve allow the patient to maximize support. The athletic shoe tread provides shock-absorbing technology and stable footing during gait. Malleable arms allow the uprights to be customized to the patient for an intimate fit and improved comfort. Available in tall and short configurations, the MAC is ideal for strains and sprains of the foot and ankle – stable foot fractures as well as distal tibial and fibular fractures. A toe cover accessory for added protection of the toes is available.

- Suggested HCPCS: OTS L4361 or CF L4360

Product Size	Foot Length	Mens Shoe Size	Womens Shoe Size
XS	Up to 9"	Up to 5	4 – 6
S	9" – 9.75"	5.5 – 7.5	6.5 – 9
M	9.75" – 10.75"	7.5 – 10.5	9.5 – 12
L	10.75" – 11.75"	10.5 – 13	12 – 14
XL	11.75" – 12.5"	13.5 – 16	14 +

MAC Achilles Wedge Insole



Product Name	Part Number
MAC Extra Small Tall Boot	MACAIR01-XS
MAC Small Tall Boot	MACAIR01-S
MAC Medium Tall Boot	MACAIR01-M
MAC Large Tall Boot	MACAIR01-L
MAC Extra Large Tall Boot	MACAIR01-XL
MAC Extra Small Short Boot	MACAIR01-S-S
MAC Small Short Boot	MACAIR01-S-S
MAC Medium Short Boot	MACAIR01-S-M
MAC Large Short Boot	MACAIR01-S-L
MAC Extra Large Short Boot	MACAIR01-S-XL
Soft Toe Cover Accessory	MACTOE01
MAC Achilles Wedge Insole	MACAWI01
Terry Cloth Liner for MAC Boot	MACLINER01



Royal Boot® Air

The Royal Boot Air is a full-shell walker boot. Features include a pneumatic bladder system to control swelling and ensure a customized fit, and a rocker-bottom to simulate natural gait.

- Suggested HCPCS: OTS L4361 or CF L4360

Product Size	Foot Length	Mens Shoe Size	Womens Shoe Size
XS	Up to 9"	Up to 5	4 – 6
S	9" – 9.75"	5.5 – 7.5	6.5 – 9
M	9.75" – 10.75"	7.5 – 10.5	9.5 – 12
L	10.75" – 11.75"	10.5 – 13	12 – 14
XL	11.75" – 12.5"	13.5 – 16	14 +

Product Name	Part Number
Royal Boot Air Extra Small Short	SRB001XS
Royal Boot Air Small Short	SRB001S
Royal Boot Air Medium Short	SRB001M
Royal Boot Air L Short	SRB001L
Royal Boot Air XL Short	SRB001XL
Royal Boot Air XS	RB001XS
Royal Boot Air Small	RB001S
Royal Boot Air Medium	RB001M
Royal Boot Air L	RB001L
Royal Boot Air XL	RB001XL
Soft Toe Cover Accessory	MACTOE01



ManaEZ® Boot Air

A walking boot that helps stabilize the foot and ankle after an injury. Features include a pneumatic air bladder system to help control swelling and ensure a snug fit, and a rocker-bottom to simulate natural gait.

- Pneumatic pump with tall and short versions
- Suggested HCPCS: OTS L4361 or CF L4360

Product Size	Foot Length	Mens Shoe Size	Womens Shoe Size
XS	Up to 9"	Up to 5	4 – 6
S	9" – 9.75"	5.5 – 7.5	6.5 – 9
M	9.75" – 10.75"	7.5 – 10.5	9.5 – 12
L	10.75" – 11.75"	10.5 – 13	12 – 14
XL	11.75" – 12.5"	13.5 – 16	14 +

Product Name	Part Number
ManaEZ Boot Air Extra Small	EZBA01XS
ManaEZ Boot Air Small	EZBA01S
ManaEZ Boot Air Medium	EZBA01M
ManaEZ Boot Air Large	EZBA01L
ManaEZ Boot Air XLarge	EZBA01XL
ManaEZ Boot Air Extra Small Short	SEZBA01XS
ManaEZ Boot Air Small Short	SEZBA01S
ManaEZ Boot Air Medium Short	SEZBA01M
ManaEZ Boot Air L Short	SEZBA01L
ManaEZ Boot Air XL Short	SEZBA01XL
Soft Toe Cover Accessory	MACTOE01



ManaEZ® Boot Pediatric

ManaEZ® Boot Pediatric

The ManaEZ Boot Pediatric walking boot for children that stabilizes the foot and ankle after an injury. Features include color-coding for easy size selection, anatomically shaped plastic uprights to help facilitate a secure fit for children, universal left/right, and easy to apply and adjust with a simple strapping and liner system. Personalize your boot with the included sticker pack!

- Suggested HCPCS: L4387

Product Size	Child's Foot Length	Shoe Size
S	Up to 6"	Up to 8.5
M	6" - 7"	8.5 - 11.5
L	7" - 8"	11.5 - 1.5 (Youth)

Product Name	Part Number
ManaEZ Boot Pediatric Small (Red)	EZBPEDS
ManaEZ Boot Pediatric Medium (White)	EZBPEDM
ManaEZ Boot Pediatric Large (Blue)	EZBPEDL



ManaGO® 71

ManaGO® 71

The ManaGO 71 provides three different levels of support in one brace. The unique design allows components to be removed as the patient's condition improves. In its extended gauntlet configuration, the ManaGO 71 provides maximum protection for ankle sprains and Syndesmosis injuries. The gauntlet and hinge stirrup with footplate captures the entire ankle mechanism. Removable upper cuff and hinge stirrup with footplate can be removed as recovery improves.

Features:

- Speed lacing ankle foot brace that extends to boot height on the calf
- Interlocking plastic plates and hinge assembly
- Figure-eight strapping

Product Size	Men's Shoe Size	Women's Shoe Size	Ankle Circumference
S	6.5 - 8	7.5 - 9	12" - 13"
M	8.5 - 10	9.5 - 11	13" - 14"
L	10.5 - 12	11.5 - 13	14" - 15"
XL	12.5 - 14	13.5 - 15	15" - 16"

Product Name	Part Number
ManaGO 71 Brace Small	MG071S
ManaGO 71 Brace Medium	MG071M
ManaGO 71 Brace Large	MG071L
ManaGO 71 Brace Extra Large	MG071XL



ManaEZ® Ankle 71

ManaEZ® Ankle 71

The ManaEZ Ankle 71 provides maximum control of the foot and ankle without total immobilization. Features a rigid plastic cuff, hinge assembly, and a lacer gauntlet for a secure fit.

Features:

- Lightweight with a breathable liner, step down to a simple figure-eight ankle brace
- Suggested HCPCS: L1971

Product Name	Part Number
ManaEZ Ankle 71 Left Small	EZA71LS
ManaEZ Ankle 71 Right Small	EZA71RS
ManaEZ Ankle 71 Left Medium	EZA71LM
ManaEZ Ankle 71 Right Medium	EZA71RM
ManaEZ Ankle 71 Left Large	EZA71LL
ManaEZ Ankle 71 Right Large	EZA71RL



ManaEZ® Ankle 06

ManaEZ® Ankle 06

The Ankle 06 incorporates an anatomical hinge with the figure 8 ankle gauntlet. A unique universal left/right hinge design ensures optimum support. The footplate provides semi-rigid support from the heel to the midfoot without restricting plantar flexion or extension. Ideal for patients with ankle strains and sprains. The unique design provides support from the acute injury phase all the way to prophylactic support during activities.

- Suggested HCPCS: L1906

Product Name	Part Number
ManaEZ Ankle 06 Extra Small	EZANKLE06XS
ManaEZ Ankle 06 Small	EZANKLE06S
ManaEZ Ankle 06 Medium	EZANKLE06M
ManaEZ Ankle 06 Large	EZANKLE06L
ManaEZ Ankle 06 Extra Large	EZANKLE06XL
ManaEZ Ankle 06 Extra Extra Large	EZANKLE06XXL



ManaEZ Ankle 02

ManaEZ® Ankle 02

A low-profile multi-ligament ankle orthosis.

Features:

- Speed lacer closure for quick application
- Breathable material, latex-free, and lightweight
- Figure-eight ligamentous support

- Suggested HCPCS: L1902

Product Name	Part Number
ManaEZ Ankle 02 Extra Small	EZANKLE02XS
ManaEZ Ankle 02 Small	EZANKLE02S
ManaEZ Ankle 02 Medium	EZANKLE02M
ManaEZ Ankle 02 Large	EZANKLE02L
ManaEZ Ankle 02 Extra Large	EZANKLE02XL
ManaEZ Ankle 02 Extra Extra Large	EZANKLE02XXL



ManaEZ Ankle Stirrup

ManaEZ® Ankle Stirrup

An easy-to-apply rigid ankle orthosis. Features adjustable straps that accommodate universal sizing for the left or right foot and a low profile for minimal bulk in the shoe. Used for acute injury through rehabilitation.

- Suggested HCPCS: L4350

Product Name	Part Number
ManaEZ Ankle Stirrup (Universal)	EZAS4350



ManaEZ Post-Op
SQ Toe Shoe

ManaEZ® Post-Op Shoe

Ideal for patients who need support and protection after trauma or post-op where a walker boot is too much immobilization, and the patient cannot fit into a standard shoe. A durable tread design provides stability during gait. Easy-to-apply universal left/right design.

- Suggested HCPCS: L3260

Product Name	Part Number
ManaEZ Post-Op Shoe Small	EZPS01-S
ManaEZ Post-Op Shoe Medium	EZPS01-M
ManaEZ Post-Op Shoe Large	EZPS01-L
ManaEZ Post-Op Shoe Extra Large	EZPS01-XL



ManaEZ Night Splint
Plantar Fasciitis

ManaEZ® Night Splint Plantar Fasciitis

Ideal for patients that need to stretch the Achilles Tendon or Plantar Fascia. Common indications are Plantar Fasciitis, Tight Achilles, Achilles Tendonitis and conditions where dorsiflexing the foot is part of the rehabilitation protocol. Easy-to-apply universal left/right design available in three sizes.

- Suggested HCPCS: L3260

Product Name	Part Number
ManaEZ Night Splint Small	EZSPNT01-S
ManaEZ Night Splint Medium	EZSPNT01-M
ManaEZ Night Splint Large	EZSPNT01-L

Spine





Tailback 2-64

Tailback 2-64

The Tailback 2-64 is a universal spinal orthosis designed with our dual-pull system. It provides a rigid, comfortable, and easy-to-apply solution for every patient. Ice pack included.

Features:

- Universal TLSO with adjustable sternal plate (25"– 64")
- Moldable and rigid posterior, anterior, and sagittal panels
- Posterior panel extends from scapular spine to sacrococcygeal junction, and anterior panel extends from pubis to sternal notch
- Shoulder straps assist in restricting gross trunk motion
- Dual-pull system
- Lightweight, latex-free, and breathable material with step-down capabilities

- Suggested HCPCS: L0464

Product Name	Part Number
Tailback 2-64 Brace OTS (Size Universal)	TB064-2
Tailback 64 Brace OTS (Single-pull)	TB064
Tailback Extension Strap	TBEXT
Tailback Ice Pack	TBICE



Tailback 2-56

Tailback 2-56

The Tailback 2-56 is a universal spinal orthosis designed with our dual-pull system. It provides a rigid, comfortable, and easy-to-apply solution for every patient. Ice pack included.

Features:

- Universal TLSO (25"– 64")
- Moldable and rigid posterior, anterior, and sagittal panels
- Posterior panel extends from scapular spine to sacrococcygeal junction, soft anterior apron
- Shoulder straps assist in maintaining proper positioning
- Dual-pull system
- Lightweight, latex-free, and breathable material with step-down capabilities

- Suggested HCPCS: L0456

Product Name	Part Number
Tailback 2-56 Brace OTS (Size Universal)	TB056-2
Tailback 56 Brace OTS (Single-pull)	TB056
Tailback Extension Strap	TBEXT
Tailback Ice Pack	TBICE



Tailback 2-50

Tailback 2-50

The Tailback 2-50 is a universal spinal orthosis designed with our dual-pull system. It provides a rigid, comfortable, and easy-to-apply solution for every patient. Ice pack included.

Features:

- Moldable and rigid posterior, anterior, and sagittal panels
- Extends from T9 to sacrococcygeal junction
- Dual-pull system
- Universal sizing (25"– 64"), latex-free, lightweight, and breathable with step-down capabilities

- Suggested HCPCS:
OTS L0650 or CF L0637

Product Name	Part Number
Tailback 2-50 Brace (Size Universal)	TB050-2
Tailback 50 Brace (Single-pull)	TB050
Tailback Extension Strap	TBEXT
Tailback Ice Pack	TBICE



Tailback 2-48

Tailback 2-48

The Tailback 2-48 is a universal spinal orthosis designed with our dual-pull system. It provides a rigid, comfortable, and easy-to-apply solution for every patient. Ice pack included.

Features:

- Rigid anterior and posterior panels
- Extends from T9 to sacrococcygeal junction
- Provides sagittal control
- Dual-pull system
- Universal sizing (25"– 64"), latex-free, lightweight, and breathable with step-down capabilities

- Suggested HCPCS: OTS L0648 or CF L0631

Product Name	Part Number
Tailback 2-48 Brace (Size Universal)	TB048-2
Tailback 48 Brace (Single-pull)	TB048
Tailback Extension Strap	TBEXT
Tailback Ice Pack	TBICE



Tailback 2-42

Tailback 2-42

The Tailback 2-42 is a universal spinal orthosis designed with our dual-pull system. It provides a rigid, comfortable, and easy-to-apply solution for every patient. Ice pack included.

Features:

- Rigid anterior and posterior panel
- Extends from L5 to L1
- Dual-pull system
- Universal sizing (25"– 64"), latex-free, lightweight, and breathable with step-down capabilities

- Suggested HCPCS: OTS L0627 or CF L0642

Product Name	Part Number
Tailback 2-42 Brace OTS (Size Universal)	TB042-2
Tailback 42 Brace OTS (Single-pull)	TB042
Tailback Extension Strap	TBEXT
Tailback Ice Pack	TBICE



Tailback 27

Tailback 27

The Tailback 27 is a lumbar orthosis that is low profile with anterior and posterior panels. The dual-pull closure makes brace adjustments quickly and easily. Soft breathable material conforms to the body for an intimate fit that can be worn over or under loose-fitting clothes. Every brace comes with an ice pack to reduce swelling and pain. Ideal for strains and sprains of the lumbar spine, degenerative disc disease, bulging disc, laminectomy, and facet syndrome.

Features:

- Wrap-around design with hook and loop belt loops
- Medium-density foam body is latex-free, low-profile
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance

- Suggested HCPCS: OTS L0642 or CF L0627

Product Name	Part Number
Tailback 27	TB027-2



Halfback 50

Halfback 50

Halfback50 Lumbar-Sacral Orthosis (LSO) provides support for spinal surgeries and conservative treatments of spinal injuries. It provides a rigid, comfortable, and easy-to-apply solution for every patient.

Features:

- Moldable and rigid posterior, anterior, and sagittal panels
- Extends from T9 to sacrococcygeal junction
- Single-pull system
- Universal sizing (25"– 64"), latex-free lightweight, and enhanced breathability with step-down capabilities

- Suggested HCPCS: OTS L0650 or CF L0637

Product Name	Part Number
Halfback 50	HB050



ManaMomma® Maternity Back Brace

A universal spinal orthosis designed to provide support to patients with the pendulous abdomen. Experience peace of mind and comfort. ManaMomma is designed to support you and your baby every step of the way. Specifically designed to limit anterior pressure on the fetus, ensuring maximum safety and comfort for both mother and baby.

Features:

- Easily adjusted throughout pregnancy, accommodating your growing belly and providing consistent support.
- Made from soft, breathable materials that are gentle on the skin, ensuring you stay comfortable all day long.
- Simple to wear and adjust, making it a hassle-free addition to your daily routine.

- Suggested HCPCS:OTS L0650 or CF L0637

Product Name	Part Number
ManaMomma Maternity Back Brace	MM001



ManaEZ® Collar 80

The ManaEZ Collar 80 is a semi-rigid thermoplastic collar, ideal for strains and sprains of the cervical vertebrae. Two-piece design with thoracic extension.

Features:

- Adjustable and lockable mandibular and occipital support
- Lightweight, latex-free, washable and breathable lining

- Suggested HCPCS:OTS L1499 or CF L0180

Product Name	Part Number
ManaEZ Collar 80 (Size Universal)	EZC080



ManaEZ® Collar 74

The ManaEZ Collar 74 is a semi-rigid thermoplastic collar, ideal for strains and sprains of the cervical vertebrae. Two-piece design with thoracic extension.

Features:

- Adjustable and lockable mandibular support
- Lightweight, latex-free, washable and breathable lining

- Suggested HCPCS: OTS L1499 or CF L0174

Product Name	Part Number
ManaEZ Collar 74 (Size Universal)	EZC074



ManaEZ® Soft Collar Universal

ManaEZ Soft Collar Universal provides uniform comfort and support for mild strains and sprains of the cervical spine. Universally sized, lightweight, breathable, and latex-free with a 3.5" medium-density foam.

- Suggested HCPCS: L0120

Product Name	Part Number
ManaEZ Soft Collar (Size) Universal	EZ120

Pain and Recovery Solutions



ManaSport+

The ManaSport+ is a portable and rechargeable prescriptive LIPUS device (low-intensity pulsed ultrasound) that is ideal for treating pain. The device generates heat within the body tissues for the treatment of muscle spasms and joint contractures, increasing local circulation and pain relief.

Features:

- Comes with a carrying case, charger, two tubes of coupling gel, reusable transducer applicator, and a transducer mounting strap
- 20-minute wear time
- Compliance tracking that shows up to 3x/day use on a calendar
- Ideal for soft tissue injuries
- Rechargeable battery with a 4-hour charge time, with 10 hours of use on full battery
- Wireless for easy patient use
- Suggested HCPCS: K1004 or E1399



ManaSport+

Product Name	Part Number
ManaSport+	MSPORT02
Ultrasound Coupling Gel	MANAGELO1
Transducer Fixation Strap	MTFIXSTRAP
Transducer Holder Patch	MTHPATCH
Transducer Weighted Strap	MTWTSTRAP
Ultrasound Replacement Charger	MSUSCHG



ManaFlexx 2

ManaFlexx 2

- FDA cleared NMES (mode 1) and TENS (mode 2) device
- Comes with initial electrode pad
- Provides muscle re-education and strengthening
- Maintains or increases range of motion and reduces muscle spasms
- Immediate post-surgical stimulation of calf muscles to prevent venous thrombosis
- Additional lead wires available to extend your reach

- Suggested HCPCS: NMES E0745, TENS E0720



ManaFlexx 2 Electrode Sheet 10 pack



ManaFlexx 2 Lead Wires

Product Name	Part Number
ManaFlexx 2	MF002-RX
ManaFlexx 2 Replaceable Pad	MFPAD
ManaFlexx 2 Lead Wires	MF002LW
ManaFlexx 2 Electrode Sheet 10 pack	MFELP10



ManaCold® Platinum 2.0 Universal

ManaCold Platinum 2.0 Universal

ManaCold Platinum 2.0 Universal Combo Unit is a water-circulating therapy device used to decrease pain, swelling, and muscle spasms after surgery or injury. When the cold therapy pad is applied to the patient's injured area, it delivers continuous cold therapy through the pad. It creates an easy solution for your patients without a large mess.

Features:

- Contains a 4-quart cooling reservoir, 7' insulation tube, and a cold therapy pad
- Latex-free, compact, and easy to store
- ManaCold Platinum Universal can be used and operated by a patient or healthcare provider in a home or clinical setting
- Plugs into wall during use
- Features a leak-proof lid and our most secure connectors
- Hand control for customized prescription
- Joint specific pads

• Unit Suggested HCPCS: E0218

• Universal Pad Suggested HCPCS: E1399 or E0249



ManaCold Platinum
Universal Circulating Pad



ManaCold Platinum
Ankle Pad



ManaCold Platinum
Back Pad



ManaCold Platinum
Hip Pad



ManaCold Platinum
Shoulder Pad



Product Name	Part Number
ManaCold Platinum 2.0 Universal Combo Unit	MANACOLD02-PT
ManaCold Platinum 2.0 Universal Pad	MCOLDPAD02
Ankle Pad 2.0	MCOLDPAD02-ANKLE
Back Pad 2.0	MCOLDPAD02-BACK
Hip Pad 2.0	MCOLDPAD02-HIP
Knee Pad 2.0	MCOLDPAD02-KNEE
Shoulder Pad 2.0	MCOLDPAD02-SHLD
ManaCold Platinum Replacement Charger	MCOLDCHG

Product Name	Part Number
ManaCold Platinum Universal Combo Unit	MANACOLD01-PT
ManaCold Platinum Universal Pad	MCOLDPAD
Ankle Pad	MCOLDPAD-ANKLE
Back Pad	MCOLDPAD-BACK
Hip Pad	MCOLDPAD-HIP
Shoulder Pad	MCOLDPAD-SHLD
ManaCold Platinum Replacement Charger	MCOLDCHG

Braces with Ice and Compression



ManaEZ[®] Hip Ice

ManaEZ[®] Hip Ice

The ManaEZ Hip Ice is a hip abduction control orthosis, which allows protected motion while providing edema and pain relief with ice and compression.

Features:

- Patented ROM hinge consists of an abduction pivot hinge and a ROM flexion/extension hinge to protect the hip from every angle
- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance
- Circumferential straps on waist and thigh for easy application
- Universally sized for left or right hip; latex-free

- Suggested HCPCS: OTS L2999 or CF L1686

Product Name	Part Number
ManaEZ Hip Ice	EZHIP01B
ManaEZ Hip Ice Extension Strap	EZIWHIPEXT



ManaEZ[®] ROM Ice

ManaEZ[®] ROM Ice

The ManaEZ ROM Ice is a hinged brace that provides a protected range of motion with the added benefit of ice and compression. Ideal for patients on the move.

Features:

- ROM hinges set a protected range of motion with a safety cap
- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance
- Posterior wrap-around design for easy application
- Universally sized for left or right leg; latex-free

- Suggested HCPCS: OTS L1833 or CF L1832

Product Name	Part Number
ManaEZ ROM Ice (Size Universal)	EZR001
ManaEZ ROM Ice Extra Small Size	EZR001-XSM
ManaEZ ROM Ice Pack	EZR001-IP
ManaEZ ROM Ice Extension Strap	EZROMEXT



ManaEZ[®] Hinged Ice

ManaEZ[®] Hinged Ice

The ManaEZ Hinged Ice is a quick-fitting, universally sized, poly-centric hinge orthosis which is ideal for patients who need medial/lateral support without ROM control.

Features:

- Medium-density foam structure with hinges
- Universally sized wrap-around feature for easy on/off application
- Open popliteal
- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance

- Suggested HCPCS: L1820

Product Name	Part Number
ManaEZ Hinged Ice	EZHICE



ManaEZ[®] Ankle Ice

ManaEZ[®] Ankle Ice

The ManaEZ Ankle Ice is an ankle gauntlet that controls edema and pain with ice and compression.

Features:

- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance
- Easy to apply strap closure
- Anterior closure allows for easy application
- Universally sized for left or right ankle; latex-free

- Suggested HCPCS: L1902

Product Name	Part Number
ManaEZ Ankle Ice (Size Universal)	EZANKLE01

Braces with Ice and Compression *(continued)*



ManaEZ® Shoulder Ice

ManaEZ® Shoulder Ice

The ManaEZ Shoulder Ice is a shoulder orthosis with a medium-density foam structure, which holds the shoulder in place while providing edema and pain relief with ice and compression.

Features:

- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance
- Easy to apply circumferential thoracic and bicep straps
- Universally sized for left or right shoulder; and latex-free

- Suggested HCPCS: L3670

Product Name	Part Number
ManaEZ Shoulder Ice	EZSHLD



ManaEZ® Ice Wrist

ManaEZ® Ice Wrist

The ManaEZ Ice Wrist brace provides wrist stabilization with the additional benefit of an ice pack and compression system to reduce pain and swelling. The compression is adjustable by the user. Ideal for wrist sprains, wrist surgeries, and instability of the wrist.

Features:

- Removable and malleable palmer stay
- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance
- Universally sized for left or right

- Suggested HCPCS: L3908

Product Name	Part Number
ManaEZ Ice Wrist	EZIWRIST01

Ice Wraps



ManaEZ® Ice Wrap Hip

ManaEZ® Ice Wrap Hip

The ManaEZ Ice Wrap Hip is a simple and convenient edema and pain solution for the hip and buttock.

Features:

- Wrap-around thigh and waist straps maintain ice position and ensure an easy application
- Includes 2 flexible ice packs
- Universally sized for left or right hip
- Latex-free

Product Name	Part Number
ManaEZ Ice Wrap Hip	EZIWH01
ManaEZ Ice Wrap Hip Extension Strap	EZIWHIPEXT



ManaEZ® Ice Wrap Knee

ManaEZ® Ice Wrap Knee

The ManaEZ Ice Wrap Knee is a simple but convenient edema and pain solution. Maintains ice position during ambulation.

Features:

- Ideal for knee ligament damage, knee surgeries, and instability of the knee
- Includes a flexible ice pack for pain and swelling
- It can be used for cold and heat therapy
- Latex-free

Product Name	Part Number
ManaEZ Ice Wrap Knee	EZIWKNO1

Mobility Aids



ManaEZ® Crutch

ManaEZ® Crutch

The ManaEZ Crutches are universally sized and easy to adjust with quick-release clips. Ideal for individuals from 4'6" to 6'6", up to 300 lbs.

Features:

- Comes in pairs
- High-strength, powder-coated aluminum
- Latex-free palm and arm pads
- Rubber non-slip tips

- Suggested HCPCS: E0114

Product Name	Part Number
ManaEZ Crutch (Universal)	EZCRUTCH01



ManaEZ® Crutch Spring

ManaEZ® Crutch Spring

The ManaEZ Spring Crutches are universally sized and easy to adjust with quick-release clips. These are ideal for individuals from 4'6" to 6'6", up to 300 lbs.

Features:

- Comfortable shock-absorbing spring
- Comes in pairs
- High-strength, powder-coated aluminum
- Latex-free palm and arm pads
- Rubber non-slip tips

- Suggested HCPCS: E0117 x 2

Product Name	Part Number
ManaEZ Crutch Spring (Universal)	EZSPRING01

Maternity



**ManaMomma®
Maternity Back Brace**

ManaMomma® Maternity Back Brace

A universal spinal orthosis designed to provide support to patients with the pendulous abdomen. Experience peace of mind and comfort. ManaMomma is designed to support you and your baby every step of the way. Specifically designed to limit anterior pressure on the fetus, ensuring maximum safety and comfort for both mother and baby.

Features:

- Easily adjusted throughout pregnancy, accommodating your growing belly and providing consistent support.
- Made from soft, breathable materials that are gentle on the skin, ensuring you stay comfortable all day long.
- Simple to wear and adjust, making it a hassle-free addition to your daily routine.

- Suggested HCPCS:OTS L0650 or CF L0637

Product Name	Part Number
ManaMomma Maternity Back Brace	MM001



Three-In-One Belly Band

Three-In-One Belly Band

The Three-In-One Belly Band is an elastic abdominal support that wraps around the abdomen and fastens with hook and loop. The unique design fits snugly under the belly to lift and support the abdomen. Provides abdominal and back support post-pregnancy, as well as abdominal support after abdominal or pelvic surgery. Adjustable tension through elastic materials allows customized support of the belly and lower back.

Features:

- Breathable elastic material
- Step down capabilities
- Hook and loop closure for easy application

- Suggested HCPCS:OTS L0621

Product Name	Part Number
Three-In-One Belly Band (Size Universal)	BB0301

