



MANAMED®

Recovery activated



2026

The Company

Welcome to the new chapter of ManaMed, where pioneering orthopedic solutions meet the spirit of 'Recovery Activated', driving us to the forefront of patient care and mobility. Our comprehensive product portfolio, spanning vascular therapies, assisted mobility, and DVT prevention, is designed to streamline the rehabilitation process for those eager to reclaim their active lifestyles. Our unwavering commitment to innovation keeps us constantly exploring, learning, and pushing the boundaries of what's possible in orthopedic care.

Central to our achievements are the significant partnerships we've nurtured within the industry, enhancing the level of care and offering products that expedite the return to normalcy for patients. Highlighting this is our esteemed collaboration with the UFC Performance Institute, ensuring that UFC athletes have access to the highest quality orthopedic braces, soft goods, and essential stabilizing equipment like crutches, knee braces, and walking boots. Similarly, our dedication to our distributor network—a robust community committed to excellence—reflects our ongoing investment in manufacturing, quality control, and customer service excellence.

As ManaMed continues to expand, we remain your steadfast ally in orthopedic recovery, always ready to face new challenges with creative solutions. This journey is not just about us; it's about fostering enduring connections with our customers, who are the heart and soul of our company. Our success is deeply rooted in our ability to truly understand and satisfy our customers' needs. When you partner with ManaMed, you become an integral part of a mission dedicated to promoting healthier, happier lives.



"At ManaMed, we believe that our customers are the heart and soul of our company. Our mission is to design, develop, and distribute orthopedic solutions that keep active, high-performing individuals in the game and back on track following injury or surgery. Our extensive catalog includes vascular and electrotherapies, assisted mobility, electrical stimulation, and DVT prevention products. By partnering with orthopedists, physical therapists, surgeons, and DME companies across the US, we ensure that patients of all ages and backgrounds remain healthy and comfortable on their path to recovery. When you choose us, you're not just a customer - you're part of our mission to help people live healthier, happier lives."

Trevor Theriot
President and CEO

Online Ordering Portal

**Our one stop shop for ordering,
paying invoices, and tracking orders.**

View your order history and create
ordering templates for easy reordering.
Quickly pay invoices via ACH, credit card,
or apply your credit memos.

Contact us at manamed.com/become-a-customer
and we will get you set up.

Custom Logo Program

**Add your logo to the products you order on
our portal to create better brand visibility.**


Available for most products, contact us
at logoprogram@manamed.com for
more information.



Table of Contents

Vascular Therapy	PlasmaFlow X	5	Foot and Ankle	FlexSafe Ankle Recovery Orthosis	26
	PlasmaFlow	5		MAC Boot	26
	PlasmaFlight	5		ManaEZ Boot Air 2	26
	PlasmaWave	6		ManaGO 71	27
	ManaEZ Compression Sock	6		ManaEZ Ankle 71	27
				ManaEZ Ankle 06	27
				ManaEZ Ankle 02	27
Lymphedema Therapy	ManaFlow 2A	8		ManaEZ Ankle Stirrup	28
	ManaFlow 2B	8		ManaEZ Post-Op SQ Toe Shoe	28
	ManaFlow 51	8		ManaEZ Night Splint	28
	ManaFlow 52	8			
Circulating Cold Therapy Devices	ManaCold Platinum 2.0 Universal	10	Spine	Tailback 64	30
	Cold Compression Wrap Pro	10		Tailback 56	30
Pain and Recovery Solutions	ManaRay	11		Tailback 2-50	30
	ManaSport+	12		Tailback 48	31
	ManaFlexx 2	13		Tailback 2-42	31
				Halfback 2-50	31
				Tailback 27	32
Braces with Ice and Compression	ManaEZ Hip Ice	14		Tailback Extension Strap	32
	ManaEZ ROM Ice	14		ManaEZ Ice Long Duration Back Pad	32
	ManaEZ Hinged Ice	14		TB050 Ice Pack	32
	ManaEZ Ankle Ice	15		ManaEZ Collar 80	33
	ManaEZ Shoulder Ice	15		ManaEZ Collar 74	33
	ManaEZ Ice Wrist	15		ManaEZ Soft Collar	33
Ice Wraps	ManaEZ Ice Wrap Hip	16	Surgical Bundles and Supplies	Total Joint Bundle	35
	ManaEZ Ice Wrap Knee	16		Spine Bundle	35
	ManaEZ Ice Wrap Wrist	16		Ankle Bundle	35
				Hip Bundle	35
Shoulder	Kahuna Brace	18		Shoulder Bundle	36
	Kahuna 3960 Shoulder Sling	18		Cervical Bundle	36
	Mini Kahuna Brace	18		ManaMomma Maternity Bundle	36
	ManaEZ Sling	18		Reddy Joint Bath	37
	ManaEZ Humerus	13		ManaSlide	37
	Mana Water Sling	13			
Elbow	Cabo Elbow	19	Maternity	ManaMomma Maternity Back Brace	37
				Three-In-One Belly Band	38
Hand and Wrist	ManaEZ Thumb Spica	20	Mobility Aids	ManaEZ Crutch	38
	ManaEZ Wrist	20		ManaEZ Crutch Spring	38
			Online Access	Device Videos	39
Knee	ManaEZ ROM Elite/ManaEZ ROM Elite Short	22		Fitting Instructions	39
	ManaEZ Hinged Knee	22			
	ManaEZ Wrap 33	22			
	ManaEZ Knee Immobilizer	23			
	EZOA	23			
	GOAT	23			
	Cabo Knee	24			
	ManaEZ Tendon Strap	24			
	ManaEZ Suspension Sleeve	24			

Vascular Therapy





PlasmaFlow® X

PlasmaFlow® X

The new PlasmaFlow® X is intended to be an easy-to-use sequential compression system prescribed by a physician, for use in the home or clinical setting to help prevent the onset of DVT in patients by simulating muscle contractions. This compact device enhances blood circulation, diminishing post-operative pain and swelling, and enhancing wound healing.

Features:

- Extended battery life – up to 19 hours on a fast and reliable USB-C charge
- Dual inflation modes – Continuous and Step-Up inflation up to 55 mmHg
- Large LCD display – shows usage, pressure, battery status, and mode
- LCD light sleep mode – dims the display to avoid disturbing patients at night
- Innovative calf-capturing design – secure and comfortable fit
- Compact design – approximately 30% reduction in shelf space
- Portable and tubeless – lightweight for home or travel

- Suggested HCPC Code(s): E0676

Product Name	Part Number
PlasmaFlow X	XPF0001
Accessories	Part Number
Extension Sleeves	XPF-EXT
PlasmaFlow X Replacement Charger	XPFCHG



PlasmaFlow®

PlasmaFlow®

The PlasmaFlow is intended to be an easy-to-use sequential compression system prescribed by a physician. For use in the home or clinical setting to help prevent the onset of DVT in patients, by stimulating blood flow in the extremities by simulating muscle contractions. This device enhances blood circulation, diminishing post-operative pain and swelling, and enhancing wound healing.

Features:

- 2 Modes: Continuous inflation to 55mmHG then deflate and Step-Up inflation technology which allows the unit to increase pressure in slower 10mmHG increments to 55mmHG
- Long-lasting rechargeable battery. Up to 10 hours of use on one charge
- 2 LCD screens to monitor usage and pressure
- Portable, lightweight, and tubeless
- Suggested HCPC Code(s): E0676

Product Name	Part Number
PlasmaFlow DVT	PF0001
Accessories	Part Number
PlasmaFlow Extension Sleeves	PFEXT
PlasmaFlow Charger	PFCHG
PlasmaFlow Hospital Charger	PFHCHG
PlasmaFlow Hospital 12.5' Charger	PFCHG-12.5



PlasmaFlight®

PlasmaFlight®

The PlasmaFlight® has been approved by the FDA as an over-the-counter portable inflatable massage system that simulates the kneading and stroking of tissue with the hands by use of an inflatable cuff. This device is used for increasing blood circulation in the treated areas, as well as relief of minor muscle aches and pains.

Features:

- Pump will inflate each leg cuff to pre-set pressure of 55mmHG and releases once that pressure is reached
- Visual display shows hourly usage, mode, and pressure level
- Electronically controlled pumps compress the calf or calves to aid in blood flow in the lower extremity
- Internal rechargeable battery and cordless design allow for ease of placement and patient mobility

Product Name	Part Number
PlasmaFlight	PFLIGHT01



PlasmaWave®

PlasmaWave®

The PlasmaWave is a game-changing clinical grade pneumatic SCD device used in the home or clinical setting. It has a market-leading 25+ hour battery.

Features:

- Over 8 feet of reusable tubing
- 25+ hour battery
- Sleeve options include both reusable or disposable
- Convenient bed rack also available
- 4 Unique pressure modes which allows the unit to increase pressure in slower increments
- 1 LCD screen to monitor usage, pressure, and battery life
- Lightweight
- Suggested HCPC Code(s): E0676

Product Name	Part Number
PlasmaWave Portable DVT	PWAVE01
Compression Sleeves	
PlasmaWave Standard Calf Sleeves	M-40
PlasmaWave XL Calf Sleeves	M-42
PlasmaWave Standard Foot Garment	M-21
PlasmaWave Standard Thigh Garment	M-50
Go Green Sleeve Barrier (10 pairs per box)	M-32
M-31 Reusable Universal Calf Garment	M-31
Accessories	
PlasmaWave Bed Rack	PWRACK
PlasmaWave Bed Holder	PWBH01
PlasmaWave Tube Set	PW15
PlasmaWave Replacement Charger	PWCHG



M-31 Reusable Sleeve
• Suggested HCPC Code(s): E0673



Go Green
Compression Sleeve Barrier



"The Eco Friendly Sleeve"



ManaEZ® Compression Sock

ManaEZ Compression Sock

The ManaEZ Compression sock is a quick-fitting, sized sock that uses a gradient descent compression for up to 20-30 mmHg of compression. The sock is made with breathable nylon/spandex material that is quick-drying and has enhanced air circulation and odor fighting properties.

- Latex free
- Unique 3D stereoscopic pattern design that enhances shock absorption
- Machine washable with cold water
- Suggested HCPC Code(s): A6530

Product Name	Part Number
ManaEZ Compression Sock Medium	EZSOCK01M
ManaEZ Compression Sock Large	EZSOCK01L
ManaEZ Compression Sock Extra Large	EZSOCK01XL

Lymphedema Therapy





ManaFlow® 2A

ManaFlow® 2A

The ManaFlow 2A provides gradient sequential compression with customizable pressures in each of the 6 chambers. ManaFlow 2A can be adjusted by the physician to a pressure of 20, 30, 40, 50, 60, 70, or 80 mmHg to treat lymphedema, treatment of chronic venous insufficiency, or enhance venous return.

Features:

- Clearing therapy mode
- Fully integrated controller
- Flexible sizing — zipper closure
- USB-C charging — usable while plugged in
- 6 chamber garment — adjustable gradient
- Clear LCD — see & adjust with confidence
- Cordless power — up to 8 hours
- Easy-to-clean material

- Suggested HCPC Code(s): E0652 and E0667

Product Name	Part Number
ManaFlow 2A	MFLOW2A
Accessories	Part Number
ManaFlow 2A Replacement Charger	MFLOW2CHG
ManaFlow 2A Full Leg Sleeve	MFLOW-UFLS
ManaFlow 2A Extender Strap	MFLOW-UFLSEXT



ManaFlow® 2B

ManaFlow® 2B

The ManaFlow 2B has three different pressure options of Low (30 mmHg), Medium (40 mmHg), and High (60mmHg) which all 6 chambers are set to. ManaFlow 2B is designed to treat lymphedema, treatment of chronic venous insufficiency, or enhance venous return.

Features:

- Fully integrated controller
- Flexible sizing — zipper closure
- USB-C charging — usable while plugged in
- 6 chamber garment — adjustable gradient
- Clear LCD — see & adjust with confidence
- Cordless power — up to 8 hours
- Easy-to-clean material

- Suggested HCPC Code(s): E0651 and E0667

Product Name	Part Number
ManaFlow 2B	MFLOW2B
Accessories	Part Number
ManaFlow 2B Replacement Charger	MFLOW2CHG
ManaFlow 2B Full Leg Sleeve	MFLOW-UFLS
ManaFlow 2B Extender Strap	MFLOW-UFLSEXT



ManaFlow® 51/52

ManaFlow® 51 (Uncalibrated) and ManaFlow® 52 (Calibrated)

The ManaFlow 51 and 52 are four-chamber tubeless pneumatic compression devices, used to treat lymphedema, treatment of chronic venous insufficiency, or enhance venous return. The ManaFlow 51 provides a preset 50 mmHg in each of the four chambers, while the ManaFlow 52 allows the clinician to calibrate the pressure of each chamber.

Features:

- Both available in calf height or full leg
- Rechargeable battery takes 3 hours to charge, 8 hour use on a single charge, can be used while plugged in
- Compliance meter shows use, pressure, and chamber
- Gradient compression
- Easy-to-clean material

- Suggested HCPC Code(s): E0651, E0662, E0657 and E0669

Product Name	Part Number
ManaFlow 51 (Uncalibrated)	MFLOW51
ManaFlow 51 Full (Uncalibrated)	MFLOW51LEG
ManaFlow 52 (Calibrated)	MFLOW52
ManaFlow 52 Full (Calibrated)	MFLOW52LEG
ManaFlow 51/52 Extenders	MFLOWEXT

Pain and Recovery Solutions



 **MANAMED**

Recovery activated

Circulating Cold Therapy Devices



ManaCold® Platinum 2.0
Universal

ManaCold Platinum 2.0 Universal

ManaCold Platinum 2.0 Universal Combo Unit is a water-circulating therapy device used to decrease pain, swelling, and muscle spasms after surgery or injury. When the cold therapy pad is applied to the patient's injured area, it delivers continuous cold therapy through the pad. It creates an easy solution for your patients without a large mess.

Features:

- Contains a 4-quart cooling reservoir, 7' insulation tube, and a cold therapy pad
- Latex-free, compact, and easy to store
- Can be used and operated by a patient or healthcare provider in a home or clinical setting
- Plugs into wall during use
- Features a leak-proof lid and our most secure connectors
- Hand control for customized prescription
- Joint specific pads
- Unit Suggested HCPC Code(s): E0218
- Pad Suggested HCPC Code(s): E1399

Product Name	Part Number
ManaCold Platinum 2.0 Universal Combo Unit	MANACOLD02-PT
Wrap-On Pads	Part Number
ManaCold Platinum 2.0 Universal Pad	MCOLDPADO2
Ankle Pad 2.0	MCOLDPADO2-ANKLE
Back Pad 2.0	MCOLDPADO2-BACK
Hip Pad 2.0	MCOLDPADO2-HIP
Knee Pad 2.0	MCOLDPADO2-KNEE
Shoulder Pad 2.0	MCOLDPADO2-SHLD
Accessories	Part Number
ManaCold Platinum Replacement Charger	MCOLDCHG



ManaCold Platinum
Ankle Pad



ManaCold Platinum
Back Pad



ManaCold Platinum
Hip Pad



ManaCold Platinum
Shoulder Pad



ManaCold Platinum
Universal Pad



Circulating Cold Therapy Devices *(continued)*



Cold Compression Wrap Pro

Cold Compression Wrap Pro

Pneumatic compression pump that attaches to one of our six compression wraps for elbow, hip, knee, shoulder, back, and ankle. Pump comes with a battery, screwdriver to insert the battery, and a charger. The compression wrap attachments are sold separately and have one removable cryo-gel pack.

Features:

- 2 pneumatic ports to treat 1 or 2 areas at one time
- Adjustable pressure settings: 20, 40, 60, 80, or 100 mmHg to control edema
- LCD screen displays pressure to show inflation cycle
- Set treatment times: 10 minute intervals up to 120 minutes with default of 30 minutes
- Has a rechargeable battery for portability or use while charging
- Easy-to-clean nylon material
- Simple strapping system for donning & doffing
- Suggested HCPC Code(s): E1399

Product Name	Part Number
Cold Compression Wrap Pro Pump	CCW-PP
Wraps	Part Number
Cold Wrap Ankle	CCW-PAW
Cold Wrap Back	CCW-PBW
Cold Wrap Elbow	CCW-PEW
Cold Wrap Hip	CCW-PHW
Cold Wrap Knee	CCW-PKW
Cold Wrap Shoulder	CCW-PSW



Cold Wrap Ankle



Cold Wrap Back



Cold Wrap Elbow



Cold Wrap Hip



Cold Wrap Knee



Cold Wrap Shoulder

Pain and Recovery Solutions

ManaRay

The ManaRay device is a wearable LED-based system with adjustable hook and loop straps designed to deliver red and near-infrared light therapy. The device is powered by a built-in 3.7V lithium-ion battery and is designed for short-duration therapy sessions to provide topical heating for the purpose of elevating tissue temperature for the temporary relief of minor muscle and joint pain, stiffness, minor arthritis pain, and muscle spasm by temporarily increasing local blood circulation and/or the temporary relaxation of muscle.

It includes:

- A control unit with USB-C charging port and power indicator lights.
- A flexible LED pad with 60 + 60 LEDs emitting at 660nm, 850nm, and 940nm.

Product Name	Part Number
ManaRay Ankle	MRAY-ANK001
ManaRay Back	MRAY-BK001
ManaRay Hand	MRAY-HND001
ManaRay Knee	MRAY-KN001
ManaRay Neck	MRAY-NK001
ManaRay Shoulder	MRAY-SH001
ManaRay Universal	MRAY-UNI001



ManaRay



ManaRay Ankle



ManaRay Back



ManaRay Hand



ManaRay Knee



ManaRay Neck



ManaRay Shoulder



ManaRay Universal

ManaSport+



ManaSport+

The ManaSport+ is a portable and rechargeable prescriptive noninvasive ultrasound device that is ideal for treating pain. The device generates heat within the body tissues for the treatment of muscle spasms and joint contractures, increasing local circulation and pain relief.

Features:

- Comes with a carrying case, charger, two tubes of coupling gel, reusable transducer applicator, and a transducer mounting strap
- 20-minute wear time
- Compliance tracking that shows up to 3x/day use on a calendar
- Ideal for soft tissue injuries
- Rechargeable battery with a 4-hour charge time, with 10 hours of use on full battery
- Wireless for easy patient use
- Suggested HCPC Code(s): K1004 or E1399



Product Name	Part Number
ManaSport+	MSPORT02
Accessories	Part Number
Ultrasound Coupling Gel	MANAGEL01
Transducer Fixation Strap	MTFIXSTRAP
Transducer Holder Patch	MTHPATCH
Transducer Weighted Strap	MTWTSTRAP
Ultrasound Replacement Charger	MSUSCHG

ManaFlexx 2



ManaFlexx[®] 2

ManaFlexx[®] 2

The ManaFlexx[®] 2 is comprised of a portable electric muscle stimulator (TENS and NMES). Used for pain relief, relaxing muscle spasms, preventing muscle atrophy, increasing blood circulation to prevent venous thrombosis, and maintaining or increasing range of motion.

Features:

- FDA & PDAC cleared NMES (mode 1) and TENS (mode 2) device
- Comes with initial reusable electrode pad
- Rechargeable, portable, and sleek device
- Additional lead wires available to extend your reach with reusable electrodes
- Place it anywhere treatment is desired
- Provides muscle re-education and strengthening
- Maintains or increases range of motion and reduces muscle spasms
- Immediate post-surgical stimulation of calf muscles to prevent venous thrombosis

- Suggested HCPC Code(s): NMES E0745, TENS E0720



ManaFlexx 2 Electrode Sheet 10 pack



ManaFlexx 2 Lead Wires

Product Name	Part Number
ManaFlexx 2	MF002-RX
ManaFlexx 2 Replaceable Pad	MFPAD
ManaFlexx 2 Lead Wires	MF002LW
ManaFlexx 2 Electrode Sheet 10 pack	MFELP10



ManaEZ[®] Hip Ice

ManaEZ[®] Hip Ice

The ManaEZ Hip Ice is a hip abduction control orthosis, which allows protected motion while providing edema and pain relief with ice and compression.

Features:

- Patented ROM hinge consists of an abduction pivot hinge and a ROM flexion/extension hinge to protect the hip from every angle
- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance
- Circumferential straps on waist and thigh for easy application
- Universally sized for left or right hip; latex-free

- Suggested HCPC Code(s): OTS L2999 or CF L1686

Product Name	Part Number
ManaEZ Hip Ice	EZHIP01B
ManaEZ Hip Ice Extension Strap	EZIWHIPEXT
ManaEZ Ice Long Duration Hip Pad	EZLDIP-HIP



ManaEZ[®] ROM Ice

ManaEZ[®] ROM Ice

The ManaEZ ROM Ice is a hinged brace that provides a protected range of motion with the added benefit of ice and compression. Ideal for patients on the move.

Features:

- ROM hinges set a protected range of motion with a safety cap
- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance
- Posterior wrap-around design for easy application
- Universally sized for left or right leg; latex-free

- Suggested HCPC Code(s): OTS L1833 or CF L1832

Product Name	Part Number
ManaEZ ROM Ice (Size Universal)	EZR001
ManaEZ ROM Ice Extra Small Size	EZR001-XSM
ManaEZ ROM Ice Pack	EZR001-IP
ManaEZ ROM Ice Extension Strap	EZROMEXT
ManaEZ Ice Long Duration Knee Pad	EZLDIP-KN



ManaEZ[®] Hinged Ice

ManaEZ[®] Hinged Ice

The ManaEZ Hinged Ice is a quick-fitting, universally sized, poly-centric hinge orthosis which is ideal for patients who need medial/lateral support without ROM control.

Features:

- Medium-density foam structure with hinges
- Universally sized wrap-around feature for easy on/off application
- Open popliteal
- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance

- Suggested HCPC Code(s): L1820

Product Name	Part Number
ManaEZ Hinged Ice	EZHICE01
ManaEZ Ice Long Duration Knee Pad	EZLDIP-KN



ManaEZ® Ankle Ice

ManaEZ® Ankle Ice

The ManaEZ Ankle Ice is an ankle gauntlet that controls edema and pain with ice and compression.

Features:

- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance
- Easy to apply strap closure
- Anterior closure allows for easy application
- Universally sized for left or right ankle; latex-free

- Suggested HCPC Code(s): L1902

Product Name	Part Number
ManaEZ Ankle Ice (Size Universal)	EZANKLE01
ManaEZ Ice Long Duration Ankle Pad	EZLDIP-ANK



ManaEZ® Shoulder Ice

ManaEZ® Shoulder Ice

The ManaEZ Shoulder Ice is a shoulder orthosis with a medium-density foam structure, which holds the shoulder in place while providing edema and pain relief with ice and compression.

Features:

- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance
- Easy to apply circumferential thoracic and bicep straps
- Universally sized for left or right shoulder; and latex-free

- Suggested HCPC Code(s): L3670

Product Name	Part Number
ManaEZ Shoulder Ice	EZSHLD
ManaEZ Ice Long Duration Shoulder Pad	EZLDIP-SHLD



ManaEZ® Ice Wrist

ManaEZ® Ice Wrist

The ManaEZ Ice Wrist brace provides wrist stabilization with the additional benefit of an ice pack and compression system to reduce pain and swelling. The compression is adjustable by the user. Ideal for wrist sprains, wrist surgeries, and instability of the wrist.

Features:

- Removable and malleable palmer stay
- Pneumatic compression technology for edema and pain control – enhances cryotherapy
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance
- Universally sized for left or right

- Suggested HCPC Code(s): L3908

Product Name	Part Number
ManaEZ Ice Wrist	EZIWRIST01



ManaEZ® Ice Wrap Hip

ManaEZ® Ice Wrap Hip

The ManaEZ Ice Wrap Hip is a simple and convenient edema and pain solution for the hip and buttock.

Features:

- Wrap-around thigh and waist straps maintain ice position and ensure an easy application
- Includes 2 flexible ice packs
- Universally sized for left or right hip
- Latex-free

Product Name	Part Number
ManaEZ Ice Wrap Hip	EZIWH01
ManaEZ Ice Wrap Hip Extension Strap	EZIWHIPEXT
ManaEZ Ice Long Duration Hip Pad	EZLDIP-HIP



ManaEZ® Ice Wrap Knee

ManaEZ® Ice Wrap Knee

The ManaEZ Ice Wrap Knee is a simple but convenient edema and pain solution. Maintains ice position during ambulation.

Features:

- Ideal for knee ligament damage, knee surgeries, and instability of the knee
- Includes a flexible ice pack for pain and swelling
- It can be used for cold and heat therapy
- Latex-free

Product Name	Part Number
ManaEZ Ice Wrap Knee	EZIWKNO1



ManaEZ® Ice Wrap Wrist

ManaEZ® Ice Wrap Wrist

The ManaEZ Ice Wrap Wrist is a simple but convenient edema and pain solution while on-the-go.

Features:

- Ideal for wrist injuries and surgeries
- Includes a flexible ice pack and a pneumatic compression bulb to enhance cryotherapy
- Universally sized for left or right wrist

Product Name	Part Number
ManaEZ Ice Wrap Wrist	EZICEWW01

Upper Extremity





Kahuna Brace

Kahuna® Brace

The Kahuna Brace is a universally-sized abduction sling and is applied to the left or right side. Features:

- Multiple positioning options; the Kahuna malleable aluminum waistband with arm support allows for unlimited arm positioning including airplane splint
- Ideal for shoulder injuries and their positioning needs
- Universal sizing reduces your inventory
- Foam ball for enhanced circulation
- Step-down abduction pillow available for purchase
- Suggested HCPC Code(s): OTS L3999 or CF L3960

Product Name	Part Number
Kahuna Brace	KB001
Replacement Kahuna Pillow	KB001-P
Replacement Kangaroo Pillow	KSS001-P



Kangaroo 3960 Shoulder Sling

Kangaroo 3960 Shoulder Sling

The Kangaroo 3960 Sling is universal and can be applied to the right or left side. The breathable material allows for comfort and optimal patient experience. Kangaroo Sling comes with patent pending abduction pillow as well as the Kahuna L-bar system. Abduction pillow allows user to carry water, coffee, and a mobile phone.

- Features:
- Ability to step-down as the patient recovers
 - Multiple positioning options with the Kahuna L-B system and the Kangaroo abduction pillow
 - Ideal for shoulder injuries and surgeries
 - Foam ball for enhanced circulation
 - Suggested HCPC Code(s): OTS L3999 or CF L3960

Product Name	Part Number
Kangaroo 3960 Shoulder Sling	KKB3960
Replacement Kahuna Pillow	KB001-P
Replacement Kangaroo Pillow	KSS001-P



Mini Kahuna® Brace

Mini Kahuna® Brace

The Mini Kahuna Brace is a universally-sized abduction sling and is applied to the left or right side. Features:

- The Kahuna Pillow included is a standard black abduction pillow
- Pillow conforms to the body and is secured with a waist strap
- Ideal for shoulder injuries and their positioning needs
- Universal sizing reduces your inventory
- Foam ball for enhanced circulation
- Suggested HCPC Code(s): OTS L3999 or CF L3670

Product Name	Part Number
Mini Kahuna Brace	MKB001
Replacement Kahuna Pillow	KB001-P
Replacement Kangaroo Pillow	KSS001-P



ManaEZ® Sling

ManaEZ® Sling

The ManaEZ Sling is a shoulder immobilizer and includes a waist strap. Features:

- Soft low-profile cuff and padded shoulder strap for optimal comfort
- Waist strap ensures proper arm position
- Is used for mild sprains, strains, stable humeral fractures, and elbow injuries
- Available in sized or universal
- Suggested HCPC Code(s): L3670

Product Name	Part Number
ManaEZ Sling Small (12.4" x 6.3")	EZSL01S
ManaEZ Sling Medium (14" x 7")	EZSL01M
ManaEZ Sling Large (15.5" x 7.9")	EZSL01L
ManaEZ Sling Extra Large (17" x 8.6")	EZSL01XL
ManaEZ Sling Extra Extra Large (18.5" x 9.5")	EZSL01XXL
ManaEZ Sling Universal	EZSLUNI



ManaEZ® Humerus

ManaEZ® Humerus Universal

The ManaEZ Humerus Universal is used to stabilize humeral fractures while allowing protected motion of the arm. Features include a shoulder cap design, customizable lightweight plastic with breathable foam, and a universal left or right arm with adjustable straps that fit most arms.

Features:

- Shoulder cap design
- Lightweight and breathable
- Mold-able and trim-able
- Universal left or right arm
- Adjustable straps fit any size
- Closed-cell foam for easy cleaning

- Suggested HCPC Code(s): L3980

Product Name	Part Number
ManaEZ Humerus Universal	EZH3980



ManaWater Sling

ManaWater Sling

The ManaWater Sling is a waterproof arm sling that is used for maintaining arm position while bathing.

Features:

- One-handed adjustment
- Lightweight and waterproof
- Adjustable wrist strap fits any size
- Can be worn on both right or left arms

- Suggested HCPC Code(s): L9900

Product Name	Part Number
ManaWater Sling	MWSLING01



Cabo® Elbow

Cabo® Elbow

This universal-fit range-of-motion elbow brace is a great alternative to casting or repetitive splinting. It treats any injury in which a limited range of elbow motion is desired.

Features:

- ROM hinge features flexion setting from 0-120 degrees and extension settings from 0-90 degrees
- Telescoping length
- Drop lock with three preset locking positions (0, 15, and 30 degrees)
- Lightweight design for comfort and ease of use
- Optional hand grip

- Suggested HCPC Code(s): OTS L3761 or CF L3760

Product Name	Part Number
Cabo Elbow Universal Left	CE001L
Cabo Elbow Universal Right	CE001R

Hand and Wrist



ManaEZ[®] Thumb Spica

ManaEZ[®] Thumb Spica

The ManaEZ Thumb Spica is a universal orthosis that includes a double speed-lacer closure.

Features:

- Low profile, high-quality, breathable material
- Three malleable thumb stays
- Palmer, lateral, and radial stays
- Suggested HCPC Code(s): OTS L3809 or CF L3807
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
ManaEZ Thumb Spica Extra Small - Left	EZTS3809L-XS
ManaEZ Thumb Spica Extra Small - Right	EZTS3809R-XS
ManaEZ Thumb Spica Universal - Left	EZTS3809L
ManaEZ Thumb Spica Universal - Right	EZTS3809R
ManaEZ Thumb Spica Extra Large - Left	EZTS3809L-XL
ManaEZ Thumb Spica Extra Large - Right	EZTS3809R-XL



ManaEZ[®] Wrist

ManaEZ[®] Wrist

The ManaEZ Wrist is a high quality, universal wrist orthosis. Single speed-lacer closure makes it easier for the user to self apply the brace.

Features:

- Low profile, high-quality, breathable material
- Palmer stay and Dorsal stay
- Suggested HCPC Code(s): L3908
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
ManaEZ Wrist Extra Small - Left	EZW3908L-XS
ManaEZ Wrist Extra Small - Right	EZW3908R-XS
ManaEZ Wrist Universal - Left	EZW3908L
ManaEZ Wrist Universal - Right	EZW3908R
ManaEZ Wrist Extra Large - Left	EZW3908L-XL
ManaEZ Wrist Extra Large - Right	EZW3908R-XL

Knee





ManaEZ[®] ROM Elite Short

ManaEZ[®] ROM Elite

ManaEZ[®] ROM Elite/ManaEZ ROM Elite Short

Features:

- ManaEZ ROM Elite has telescoping and malleable hinge bars extending from 16.5"-27" (42cm-68.5 cm)
- ManaEZ ROM Elite Short version offers a fixed length of 18"
- Hinge range of motion adjustment: Flexion (0° to 120°) and extension (0° to 90°)
- Drop lock offers the ability to control and lock the hinge motion to 0°, 15°, or 30°
- Lightweight, latex-free, and breathable material
- Easy application, universal size with quick-release buckles, and fits either left or right leg

- Suggested HCPC Code(s): OTS L1833 or CF L1832

Product Name	Part Number
ManaEZ ROM Elite (Size Universal)	EZRE001
ManaEZ ROM Elite Short (Size Universal)	EZRS1833



ManaEZ[®] Hinged Knee

ManaEZ[®] Hinged Knee

The ManaEZ Hinged Knee is a convenient wrap-on hinged sleeve. The open popliteal feature allows for a comfortable fit.

Features:

- Internal hinges enhance medial and lateral support
- Condyle pads for improved fit and added comfort
- High-quality Neoprene
- Open Popliteal

- Suggested HCPC Code(s): OTS L2999 or CF L1820

- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
ManaEZ Hinged Knee Small	EZHK1820S
ManaEZ Hinged Knee Medium	EZHK1820M
ManaEZ Hinged Knee Large	EZHK1820L
ManaEZ Hinged Knee Extra Large	EZHK1820XL
ManaEZ Hinged Knee Extra Extra Large	EZHK1820XXL
ManaEZ Hinged Knee Extra Extra Extra Large	EZHK1820XXXL



ManaEZ[®] Wrap 33

ManaEZ[®] Wrap 33

The ManaEZ Wrap 33 has unique anterior and posterior body closures, which allow the user to maintain proper hinge alignment while fitting.

Features:

- Universal sizing for left or right knee, reducing inventory
- Range of motion hinges for controlled motion
- Medial and lateral hinge provides excellent support
- Breathable neoprene body provides warmth and compression

- Suggested HCPC Code(s): OTS L1833 or CF L1832

Product Name	Part Number
ManaEZ Wrap 33 Universal	EZW33



Single Panel

3 Panel

ManaEZ®
Knee Immobilizer

ManaEZ® Knee Immobilizer

The soft, breathable design accommodates most leg circumferences, both styles can be adjusted and trimmed to provide a customized fit. The combination of the posterior stay and two adjustable lateral stays maintains the optimum leg position. The popliteal pad ensures an intimate fit around the knee joint. Available in 18", 20", and 22" lengths in our single and three-panel design. Ideal for strains and sprains of the knee, stable fractures of the tibia, and patellar tendon repairs.

- Suggested HCPC Code(s): L1830

Product Name	Part Number
ManaEZ Single Panel Knee Immobilizer 18"	EZIMM01-18
ManaEZ Single Panel Knee Immobilizer 20"	EZIMM01-20
ManaEZ Single Panel Knee Immobilizer 24"	EZIMM01-24
ManaEZ 3 Panel Knee Immobilizer 16"	EZ3IMM01-16
ManaEZ 3 Panel Knee Immobilizer 18"	EZ3PIMM01-18
ManaEZ 3 Panel Knee Immobilizer 20"	EZ3PIMM01-20
ManaEZ 3 Panel Knee Immobilizer 24"	EZ3PIMM01-24



EZOA

EZOA

The EZOA is a lightweight, low-profile, single-upright OA orthosis, ideal for moderate to severe OA.

- Adjustable varus/valgus for customized pain relief
- Butterfly adjuster allows for easy fitting without tools
- Non-slip liner designed for all-day support
- Universally sized but left or right specific
- Malleable cuffs for optimum support

- Suggested HCPC Code(s): OTS L1851 or CF L1843

Product Name	Part Number
EZOA Left Knee	OFF001L
EZOA Right Knee	OFF001R



GOAT

GOAT

The GOAT is a universally sized, low-profile, lightweight option for patients with mild to moderate knee OA.

- Butterfly varus/valgus adjustment, no tools needed
- One brace does medial and lateral offloading
- Right and left specific
- Breathable and low-profile neoprene

- Suggested HCPC Code(s): OTS L1851 or CF L1843

Product Name	Part Number
GOAT Left	GOATL
GOAT Right	GOATR



Cabo Knee®

Cabo® Knee

The Cabo Knee is a dual-upright knee brace, used for support and protection during high-demand activities.

- Lightweight aluminum alloy
- Quick-fitting buckles
- Comfortable non-slip liner
- Suggested HCPC Code(s): OTS L1852 or CF L1845
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
Cabo Knee Right Knee Brace Small	TRBL45S1
Cabo Knee Left Knee Brace Small	TRBL45S2
Cabo Knee Right Knee Brace Medium	TRBL45M1
Cabo Knee Left Knee Brace Medium	TRBL45M2
Cabo Knee Right Knee Brace Large	TRBL45L1
Cabo Knee Left Knee Brace Large	TRBL45L2
Cabo Knee Right Knee Brace Extra Large	TRBL45XL1
Cabo Knee Left Knee Brace Extra Large	TRBL45XL2
Cabo Knee Right Knee Brace XXLarge	TRBL45XXL1
Cabo Knee Left Knee Brace XXLarge	TRBL45XXL2
Accessories	Part Number
Cabo Knee Condyle Replacement Pads	TRBL-CONDYLEPAD
Strap Repair Bundle Small, Right	TRBL-STRPKITSM-R
Strap Repair Bundle Small, Left	TRBL-STRPKITSM-L
Strap Repair Bundle Medium, Right	TRBL-STRPKITM-R
Strap Repair Bundle Medium, Left	TRBL-STRPKITM-L
Strap Repair Bundle Large, Right	TRBL-STRPKITLG-R
Strap Repair Bundle Large, Left	TRBL-STRPKITLG-L
Strap Repair Bundle Extra Large, Right	TRBL-STRPKITSXL-R
Strap Repair Bundle Extra Large, Left	TRBL-STRPKITSXL-L
Strap Repair Bundle Extra Extra Large, Right	TRBL-STRPKITSXXL-R
Strap Repair Bundle Extra Extra Large, Left	TRBL-STRPKITSXXL-L



ManaEZ® Tendon Strap

ManaEZ® Tendon Strap Universal

The ManaEZ Tendon Strap Universal is used to reduce tension on the patellar tendon.

- Low profile
- Neoprene strap
- Gel pad for compression
- Suggested HCPC Code(s): L1825

Product Name	Part Number
ManaEZ Tendon Strap Universal	EZTS01



ManaEZ®
Suspension Sleeve

ManaEZ® Suspension Sleeve

Designed to keep knee braces securely in place, the Suspension Sleeve delivers reliable support and all-day comfort. Durable, flexible materials help prevent brace migration during movement while maintaining a lightweight, breathable fit. Ideal for use with knee braces and orthotics during recovery or long-term condition management.

- Features:
- Enhances stability and comfort, helps prevent knee brace migration
 - Lightweight, breathable design

Product Name	Part Number
ManaEZ Suspension Sleeve	EZSUS01

Foot and Ankle



 **MANAMED**

 **MANAMED**

Recovery activated



FlexSafe
Ankle Recovery
Orthosis

FlexSafe Ankle Recovery Orthosis

The FlexSafe Ankle Recovery Orthosis is a customizable walker boot for the foot and ankle. When applied, the foot and ankle are protected by the posterior and anterior shell. The pneumatic gauntlet system provides an intimate fit as well as reduces edema. The unique design of the soft shock absorption tread system and the unique PAFA Technology (protected ankle flex/assist) ensures comfort and a normal gait. The PAFA Technology allows protected dorsiflexion and assists with toe off.

- Suggested HCPC Code(s): OTS L4361 or CF L4360
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
FlexSafe ARO Extra Small Boot	EZARO-01XS
FlexSafe ARO Small Boot	EZARO-01S
FlexSafe ARO Medium Boot	EZARO-01M
FlexSafe ARO Large Boot	EZARO-01L
FlexSafe ARO Extra Large Boot	EZARO-01XL
Soft Toe Cover Accessory	MACTOE01



MAC Boot Tall
MAC Boot Short

MAC - ManaMed Air Customizable Walking Boot

Built-in bladder system and release valve allow the patient to maximize support. The athletic shoe tread provides shock-absorbing technology and stable footing during gait. Malleable arms allow the uprights to be customized to the patient for an intimate fit and improved comfort. Available in tall and short configurations, the MAC is ideal for strains and sprains of the foot and ankle – stable foot fractures as well as distal tibial and fibular fractures. A toe cover accessory for added protection of the toes is available.

- Suggested HCPC Code(s): OTS L4361 or CF L4360
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
MAC Extra Small Tall Boot	MACAI01-XS
MAC Small Tall Boot	MACAI01-S
MAC Medium Tall Boot	MACAI01-M
MAC Large Tall Boot	MACAI01-L
MAC Extra Large Tall Boot	MACAI01-XL
MAC Extra Small Short Boot	MACAI01-S-XS
MAC Small Short Boot	MACAI01-S-S
MAC Medium Short Boot	MACAI01-S-M
MAC Large Short Boot	MACAI01-S-L
MAC Extra Large Short Boot	MACAI01-S-XL
Soft Toe Cover Accessory	MACTOE01
MAC Achilles Wedge Insole	MACAWI01

MAC Achilles Wedge Insole



ManaEZ® Boot Air 2
ManaEZ® Boot Air 2
Short

ManaEZ® Boot Air 2

A walking boot that helps stabilize the foot and ankle after an injury. Features include a pneumatic air bladder system to help control swelling and ensure a snug fit, and a rocker-bottom to simulate natural gait.

- Pneumatic pump with tall and short versions
- Suggested HCPC Code(s): OTS L4361 or CF L4360
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
ManaEZ Boot Air 2 Extra Small	EZBA02XS
ManaEZ Boot Air 2 Small	EZBA02S
ManaEZ Boot Air 2 Medium	EZBA02M
ManaEZ Boot Air 2 Large	EZBA02L
ManaEZ Boot Air 2 Extra Large	EZBA02XL
ManaEZ Boot Air 2 Short Extra Small	SEZBA02XS
ManaEZ Boot Air 2 Short Small	SEZBA02S
ManaEZ Boot Air 2 Short Medium	SEZBA02M
ManaEZ Boot Air 2 Short Large	SEZBA02L
ManaEZ Boot Air 2 Short XLarge	SEZBA02XL
Soft Toe Cover Accessory	MACTOE01



ManaGO® 71

ManaGO® 71

The ManaGO 71 provides three different levels of support in one brace. The unique design allows components to be removed as the patient's condition improves. In its extended gauntlet configuration, the ManaGO 71 provides maximum protection for ankle sprains and Syndesmosis injuries. The gauntlet and hinge stirrup with footplate captures the entire ankle mechanism. Removable upper cuff and hinge stirrup with footplate can be removed as recovery improves.

Features:

- Speed lacing ankle foot brace that extends to boot height on the calf
- Interlocking plastic plates and hinge assembly
- Figure-eight strapping

- Suggested HCPC Code(s): L1971
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
ManaGO 71 Brace Extra Small	MG071XS
ManaGO 71 Brace Small	MG071S
ManaGO 71 Brace Medium	MG071M
ManaGO 71 Brace Large	MG071L
ManaGO 71 Brace Extra Large	MG071XL



ManaEZ® Ankle 71

ManaEZ® Ankle 71

The ManaEZ Ankle 71 provides maximum control of the foot and ankle without total immobilization. Features a rigid plastic cuff, hinge assembly, and a lacer gauntlet for a secure fit.

Features:

- Lightweight with a breathable liner, step down to a simple figure-eight ankle brace

- Suggested HCPC Code(s): L1971
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
ManaEZ Ankle 71 Left Small	EZA71LS
ManaEZ Ankle 71 Right Small	EZA71RS
ManaEZ Ankle 71 Left Medium	EZA71LM
ManaEZ Ankle 71 Right Medium	EZA71RM
ManaEZ Ankle 71 Left Large	EZA71LL
ManaEZ Ankle 71 Right Large	EZA71RL



ManaEZ® Ankle 06

ManaEZ® Ankle 06

The Ankle 06 incorporates an anatomical hinge with the figure 8 ankle gauntlet. A unique universal left/right hinge design ensures optimum support. The footplate provides semi-rigid support from the heel to the midfoot without restricting plantar flexion or extension. Ideal for patients with ankle strains and sprains. The unique design provides support from the acute injury phase all the way to prophylactic support during activities.

- Suggested HCPC Code(s): L1906
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
ManaEZ Ankle 06 Extra Small	EZANKLE06XS
ManaEZ Ankle 06 Small	EZANKLE06S
ManaEZ Ankle 06 Medium	EZANKLE06M
ManaEZ Ankle 06 Large	EZANKLE06L
ManaEZ Ankle 06 Extra Large	EZANKLE06XL



ManaEZ® Ankle 02

ManaEZ® Ankle 02

A low-profile multi-ligament ankle orthosis.

Features:

- Speed lacer closure for quick application
- Breathable material, latex-free, and lightweight
- Figure-eight ligamentous support

- Suggested HCPC Code(s): L1902
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
ManaEZ Ankle 02 Extra Small	EZANKLE02XS
ManaEZ Ankle 02 Small	EZANKLE02S
ManaEZ Ankle 02 Medium	EZANKLE02M
ManaEZ Ankle 02 Large	EZANKLE02L
ManaEZ Ankle 02 Extra Large	EZANKLE02XL



ManaEZ® Ankle Stirrup

ManaEZ® Ankle Stirrup

An easy-to-apply rigid ankle orthosis. Features adjustable straps that accommodate universal sizing for the left or right foot and a low profile for minimal bulk in the shoe. Used for acute injury through rehabilitation.

- Suggested HCPC Code(s): L4350

Product Name	Part Number
ManaEZ Ankle Stirrup (Universal)	EZAS4350

ManaEZ® Post-Op Shoe

Ideal for patients who need support and protection after trauma or post-op where a walker boot is too much immobilization, and the patient cannot fit into a standard shoe. A durable tread design provides stability during gait. Easy-to-apply universal left/right design.

- Suggested HCPC Code(s): L3260
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
ManaEZ Post-Op Shoe Small	EZPS01-S
ManaEZ Post-Op Shoe Medium	EZPS01-M
ManaEZ Post-Op Shoe Large	EZPS01-L
ManaEZ Post-Op Shoe Extra Large	EZPS01-XL



ManaEZ® Post-Op SQ Toe Shoe

ManaEZ® Night Splint Plantar Fasciitis

Ideal for patients that need to stretch the Achilles Tendon or Plantar Fascia. Common indications are Plantar Fasciitis, tight Achilles, Achilles Tendonitis and conditions where dorsiflexing the foot is part of the rehabilitation protocol. Easy-to-apply universal left/right design available in three sizes.


- Suggested HCPC Code(s): L3260
- For sizing information, scan the QR code in the back of the catalog.

Product Name	Part Number
ManaEZ Night Splint Small	EZSPNT01-S
ManaEZ Night Splint Medium	EZSPNT01-M
ManaEZ Night Splint Large	EZSPNT01-L



ManaEZ® Night Splint Plantar Fasciitis

Spine





Tailback 64

Tailback 64

The Tailback 64 is a universal spinal orthosis that provides a rigid, comfortable, and easy-to-apply solution for every patient. Standard ice pack included.

Features:

- Universal single or double-pull TLSO with adjustable sternal plate (25"– 64")
- Moldable and rigid posterior, anterior, and sagittal panels
- Posterior panel extends from scapular spine to sacrococcygeal junction, and anterior panel extends from pubis to sternal notch
- Shoulder straps assist in restricting gross trunk motion
- Lightweight, latex-free, and breathable material with step-down capabilities
- Suggested HCPC Code(s): L0464

Product Name	Part Number
Tailback 2-64 Brace OTS (Double-pull)	TB064-2
Tailback 64 Brace OTS (Single-pull)	TB064
Tailback 2-64 Brace OTS (Double-pull)	TB064-2
Tailback Extension Strap	TBEXT
Tailback Ice Pack	TB050-P
TB050 Large Size Ice Pack	TB050-P-LG
ManaEZ Ice Long Duration Back Pad	EZLDIP-BK



Tailback 56

Tailback 56

The Tailback 56 is a universal spinal orthosis that provides a rigid, comfortable, and easy-to-apply solution for every patient. Standard ice pack included.

Features:

- Universal single or double-pull TLSO (25"– 64")
- Moldable and rigid posterior, anterior, and sagittal panels
- Posterior panel extends from scapular spine to sacrococcygeal junction, soft anterior apron
- Shoulder straps assist in maintaining proper positioning
- Lightweight, latex-free, and breathable material with step-down capabilities
- Suggested HCPC Code(s): L0456

Product Name	Part Number
Tailback 2-56 Brace (Double-pull)	TB056-2
Tailback 56 Brace OTS (Single-pull)	TB056
Tailback Extension Strap	TBEXT
Tailback Ice Pack	TB050-P
TB050 Large Size Ice Pack	TB050-P-LG
ManaEZ Ice Long Duration Back Pad	EZLDIP-BK



Tailback 2-50

Tailback 2-50

The Tailback 2-50 is a universal spinal orthosis designed with our dual-pull system. It provides a rigid, comfortable, and easy-to-apply solution for every patient. Standard ice pack included.

Features:

- Moldable and rigid posterior, anterior, and sagittal panels
- Extends from T9 to sacrococcygeal junction
- Dual-pull system
- Universal sizing (25"– 64"), latex-free, lightweight, and breathable with step-down capabilities

- Suggested HCPC Code(s):
OTS L0650 or CF L0637

Product Name	Part Number
Tailback 2-50 Brace (Size Universal)	TB050-2
Tailback 50 Brace (Single-pull)	TB050
Tailback Extension Strap	TBEXT
Tailback Ice Pack	TB050-P
TB050 Large Size Ice Pack	TB050-P-LG
ManaEZ Ice Long Duration Back Pad	EZLDIP-BK



Tailback 48

Tailback 48

The Tailback 48 is a universal spinal orthosis that provides a rigid, comfortable, and easy-to-apply solution for every patient. Standard ice pack included.

Features:

- Rigid anterior and posterior panels
- Extends from T9 to sacrococcygeal junction
- Provides sagittal control
- Single-pull system
- Universal sizing (25"– 64")
- Suggested HCPC Code(s): OTS L0648 or CF L0631

Product Name	Part Number
Tailback 48 Brace (Single-pull)	TB048
Tailback Extension Strap	TBEXT
Tailback Ice Pack	TB050-P
TB050 Large Size Ice Pack	TB050-P-LG
ManaEZ Ice Long Duration Back Pad	EZLDIP-BK



Tailback 2-42

Tailback 2-42

The Tailback 2-42 is a universal spinal orthosis designed with our dual-pull system. It provides a rigid, comfortable, and easy-to-apply solution for every patient. Standard ice pack included.

Features:

- Rigid anterior and posterior panel
- Extends from L5 to L1
- Dual-pull system
- Universal sizing (25"– 64"), latex-free, lightweight, and breathable with step-down capabilities
- Suggested HCPC Code(s): OTS L0627 or CF L0642

Product Name	Part Number
Tailback 2-42 Brace OTS (Size Universal)	TB042-2
Tailback Extension Strap	TBEXT
Tailback Ice Pack	TB050-P
ManaEZ Ice Long Duration Back Pad	EZLDIP-BK



Halfback 2-50

Halfback 2-50

Halfback 2-50 Lumbar-Sacral Orthosis (LSO) provides support for spinal surgeries and conservative treatments of spinal injuries. It provides a rigid, comfortable, and easy-to-apply solution for every patient. Standard ice pack included.

Features:

- Moldable and rigid posterior, anterior, and sagittal panels
- Extends from T9 to sacrococcygeal junction
- Dual-pull system
- Universal sizing (25"– 64")
- Suggested HCPC Code(s): OTS L0650 or CF L0637

Product Name	Part Number
Halfback 2-50	HB050-2
Tailback Ice Pack	TB050-P
TB050 Large Size Ice Pack	TB050-P-LG
ManaEZ Ice Long Duration Back Pad	EZLDIP-BK



Tailback 27

Tailback 27

The Tailback 27 is a lumbar orthosis that is low profile with anterior and posterior panels. The dual-pull closure makes brace adjustments quickly and easily. Soft breathable material conforms to the body for an intimate fit that can be worn over or under loose-fitting clothes. Every brace comes with a standard ice pack to reduce swelling and pain. Ideal for strains and sprains of the lumbar spine, degenerative disc disease, bulging disc, laminectomy, and facet syndrome.

Features:

- Wrap-around design with hook and loop belt loops
- Medium-density foam body is latex-free, low-profile
- Includes 2 flexible ice packs for up to two hours of cold therapy, helps ensure successful patient compliance

- Suggested HCPC Code(s): OTS L0642 or CF L0627

Product Name	Part Number
Tailback 27	TB027-2
ManaEZ Ice Long Duration Back Pad	EZLDIP-BK

Spine Accessories



Tailback Extension Strap

Tailback Extension Strap

The Tailback Extension Strap is a comfortable and easy-to-apply solution for a larger patient.

Features:

- Extends our Tailback braces an additional 12"

Product Name	Part Number
Tailback Extension Strap	TBEXT



ManaEZ[®] Ice Long Duration Back Pad

ManaEZ[®] Ice Long Duration Back Pad

Features:

- Long duration ice pack lasts 4-5x longer than a standard ice pack

Product Name	Part Number
ManaEZ Ice Long Duration Back Pad	EZLDIP-BK



TB050 Ice Pack

TB050 Large Ice Pack

Features:

- Larger than a standard ice pack

Product Name	Part Number
TB050 Large Size Ice Pack	TB050-P-LG



ManaEZ® Collar 80

ManaEZ® Collar 80

The ManaEZ Collar 80 is a semi-rigid thermoplastic collar, ideal for strains and sprains of the cervical vertebrae. Two-piece design with thoracic extension.

Features:

- Adjustable and lockable mandibular and occipital support
- Lightweight, latex-free, washable and breathable lining
- Suggested HCPC Code(s): OTS L1499 or CF L0180

Product Name	Part Number
ManaEZ Collar 80 (Size Universal)	EZC080



ManaEZ® Collar 74

ManaEZ® Collar 74

The ManaEZ Collar 74 is a semi-rigid thermoplastic collar, ideal for strains and sprains of the cervical vertebrae. Two-piece design with thoracic extension.

Features:

- Adjustable and lockable mandibular support
- Lightweight, latex-free, washable and breathable lining
- Suggested HCPC Code(s): OTS L1499 or CF L0174

Product Name	Part Number
ManaEZ Collar 74 (Size Universal)	EZC074



ManaEZ® Soft Collar

ManaEZ® Soft Collar Universal

ManaEZ Soft Collar Universal provides uniform comfort and support for mild strains and sprains of the cervical spine. Universally sized, lightweight, breathable, and latex-free with a 3.5" medium-density foam.

- Suggested HCPC Code(s): L0120

Product Name	Part Number
ManaEZ Soft Collar (Size) Universal	EZ120

Surgical Bundles and Supplies



Surgical Bundles



Total Joint Bundle

Total Joint Bundle

The Total Joint Bundle delivers a standardized post-operative solution to support pain management, mobility, and muscle re-education following joint replacement surgery.

Bundle contains:

- ManaEZ ROM Ice our cold compression device for pain, swelling, and range-of-motion control
- PlasmaFlow our Portable DVT prevention device
- ManaFlex 2 our NMES/TENS device to prevent muscular atrophy

- Suggested HCPC Code(s): L1833, E0676, E0745

Product Name	Part Number
Total Joint Bundle	TJK001



Spine Bundle

Spine Bundle

The Spine Bundle is a game changer for the DME or surgeon to provide everything your patient needs for an improved post-op outcome.

Bundle contains:

- Tailback 2-50 our LSO support
- PlasmaFlow our portable DVT prevention device
- ManaFlex 2 our NMES/TENS device to help prevent muscular atrophy

- Suggested HCPC Code(s): L0650, E0676, E0745

Product Name	Part Number
Spine Bundle	SK001



Ankle Bundle

Ankle Bundle

The Ankle Bundle is a comprehensive post-operative solution designed to support mobility, edema management, and DVT prevention following foot and ankle procedures.

Bundle contains:

- MAC Boot our pneumatic walker
- ManaEZ Ankle Ice our wrap-around ankle orthosis that has ice and compression for edema and pain relief
- PlasmaFlow our portable DVT prevention device
- ManaFlex 2 our NMES/TENS device to help prevent muscular atrophy

- Suggested HCPC Code(s): L4361, L1902, E0676, E0745

Product Name	Part Number
Ankle Bundle	AK001



Hip Bundle

Hip Bundle

The Hip Bundle delivers a standardized post-operative solution designed to support controlled motion, edema management, and patient comfort following hip surgery.

Bundle contains:

- ManaEZ ROM Hip Ice our range-of-motion cold compression wrap-around orthosis for edema and pain control
- PlasmaFlow our portable DVT prevention device
- ManaFlex 2 our NMES/TENS device to help prevent muscular atrophy

- Suggested HCPC Code(s): L1686, E0676, E0745

Product Name	Part Number
Hip Bundle	HK001



Shoulder Bundle

Shoulder Bundle

The Shoulder Bundle provides a complete post-operative solution designed to support stabilization, pain control, and recovery following shoulder surgery.

Bundle contains:

- Kahuna Brace our adjustable airplane shoulder support
- ManaEZ Shoulder Ice our wrap-around cold compression device for edema and pain control
- PlasmaFlow our portable DVT prevention device
- ManaFlexx 2 our NMES/TENS device to help prevent muscular atrophy

• Suggested HCPC Code(s): L3960, L3670, E0676, E0745

Product Name	Part Number
Shoulder Bundle	SHK001



Cervical Bundle

Cervical Bundle

The Cervical Bundle delivers a standardized post-operative solution designed to support cervical alignment, stabilization, and patient comfort following cervical spine interventions.

Bundle contains:

- ManaEZ Collar 80 our adjustable multi-post collar with mandibular and occipital support
- PlasmaFlow our portable DVT prevention device
- ManaFlexx 2 our NMES/TENS device to help prevent muscular atrophy and provide pain relief

• Suggested HCPC Code(s): L0180, E0676, E0745

Product Name	Part Number
Cervical Bundle	CK001



ManaMomma Maternity Bundle

ManaMomma Maternity Bundle

The ManaMomma Maternity Bundle provides a comprehensive support solution designed to address pre- and postpartum comfort, mobility, and recovery.

Bundle contains:

- ManaMomma Maternity Back Brace – spinal orthosis for pendulous abdomen
- Three-In-One Belly Band – elastic abdominal and pelvic support system
- ManaFlexx2 – NMES/TENS unit with a pair of lead wires and electrodes
- ManaEZ Compression Socks – gradient compression sock

• Suggested HCPC Code(s): E0745, A6530, L0650, L0621

Product Name	Part Number
ManaMomma Maternity Bundle	MMB001

Surgical Supplies



Reddy Joint Bath

Reddy Joint Bath

CDC-approved method of wound lavage with our povidone-iodine solution for PJI and SSI

Features:

- Industry-leading wound lavage
- 3-minute soak time
- Certified sterile
- Bundle includes: 5 saline solution bottles, iodine packet, bowl, and prep tray

Product Name

Reddy Joint Bath

Part Number

MMJB001B



ManaSlide

ManaSlide

ManaSlide is ideal for any patient required to be in a bed for long periods of time or need assistance in bed.

- To assist hospital staff in moving patients from bed to gurney.
- Repositioning patients that are prone to skin break down due to long duration in one position.

Product Name

ManaSlide

Part Number

MSLIDE01

Maternity



**ManaMomma®
Maternity Back Brace**

ManaMomma® Maternity Back Brace

A universal spinal orthosis designed to provide support to patients with the pendulous abdomen. Experience peace of mind and comfort. ManaMomma is designed to support you and your baby every step of the way. Specifically designed to limit anterior pressure on the fetus, ensuring maximum safety and comfort for both mother and baby.

Features:

- Easily adjusted throughout pregnancy, accommodating your growing belly and providing consistent support.
- Made from soft, breathable materials that are gentle on the skin, ensuring you stay comfortable all day long.
- Simple to wear and adjust, making it a hassle-free addition to your daily routine.

- Suggested HCPC Code(s): OTS L0650 or CF L0637

Product Name

ManaMomma Maternity Back Brace

TB050 Large Size Ice Pack

ManaEZ Ice Long Duration Back Pad

Part Number

MM001

TB050-P-LG

EZLDIP-BK



Three-In-One Belly Band

Three-In-One Belly Band

The Three-In-One Belly Band is an elastic abdominal support that wraps around the abdomen and fastens with hook and loop. The unique design fits snugly under the belly to lift and support the abdomen. Provides abdominal and back support post-pregnancy, as well as abdominal support after abdominal or pelvic surgery. Adjustable tension through elastic materials allows customized support of the belly and lower back.

Features:

- Breathable elastic material
- Step down capabilities
- Hook and loop closure for easy application

- Suggested HCPC Code(s): OTS L4467

Product Name	Part Number
Three-In-One Belly Band (Size Universal)	BB0301

Mobility Aids



ManaEZ Crutch

ManaEZ® Crutch


The ManaEZ Crutches are universally sized and easy to adjust with quick-release clips. Ideal for individuals from 4'6" to 6'6", up to 300 lbs.

Features:

- Comes in pairs
- Two high-strength, powder-coated aluminum with rubber non-slip tips
- Latex-free palm and arm pads
- Rubber non-slip tips

- Suggested HCPC Code(s): E0114

Product Name	Part Number
ManaEZ Crutch (Universal)	EZCRUTCH01



ManaEZ Crutch Spring

ManaEZ® Crutch Spring

The ManaEZ Spring Crutches are universally sized and easy to adjust with quick-release clips. These are ideal for individuals from 4'6" to 6'6", up to 300 lbs.

Features:

- Comfortable shock-absorbing spring
- Comes in pairs
- Two high-strength, powder-coated aluminum with rubber non-slip tips
- Latex-free palm and arm pads
- Rubber non-slip tips

- Suggested HCPC Code(s): E0117 x 2

Product Name	Part Number
ManaEZ Crutch Spring (Universal)	EZSPRING01

Scan a device qr code to view detailed info

Your patients simply scan the QR code with a smart phone to view a video about the ManaMed product you prescribed. The videos show how to apply and use the product, and also include tips and troubleshooting.



Scan here to view all product sizing charts

Your patients simply scan the QR code with a smartphone to see detailed sizing information on specific devices.



FlexSafe ARO



PlasmaFlow®



PlasmaWave®



ManaFlow® 51



ManaFlow® 52



ManaFlexx 2®



Mini Kahuna Brace



Bionic Sling



GOAT



PlasmaFlow® X



ManaSport +



ManaEZ® Shoulder Ice



ManaEZ® ROM Elite/
ManaEZ® ROM Elite Short



ManaEZ® Tendon Strap



ManaEZ® ROM Ice



ManaEZ® Hip Ice



ManaEZ® Collar



ManaEZ® Soft Collar



ManaEZ® Spring &
Standard Crutch



Tailback Line



Tailback 64



ManaCold® Platinum 2



ManaEZ® Ankle Stirrup



ManaEZ® Ankle 02



ManaEZ® Humerus Uni



ManaRay



Kahuna Brace



Cabo Knee®



ManaGO® 71



ManaEZ® Hinged Knee



ManaEZ® Thumb Spica



PDAC Letters

Disclaimer: The HCPCS codes listed in this catalog are provided as suggested codes based on publicly available information. We do not make claims, promises, or guarantees as to the availability of reimbursement for any product. It is within the sole discretion of the customer to determine the appropriate billing code for a product, as well as whether the use of a product complies with medical necessity and other documentation requirements of the payor. This information is not intended to be, nor should it be considered, billing or legal advice.

