

Cold Compression Wrap Pro

Cold Compression Wrap and Pneumatic Pump System

Indicated for the temporary relief of minor muscle aches and pains.



 **MANAMED**

Recovery activated

Cold Compression Wrap Pro



Indications for Use

The device is indicated for temporary increase in circulation of the treated areas in people who are in good health and simulates kneading and stroking of tissues using an inflatable wrap. The cold pack is indicated for localized therapy in situations where cold temperature therapy is necessary or desirable.

Product Description

The Cold Compression Wrap Pro delivers powerful recovery through the alliance of cold therapy and active compression. Each joint-specific garment comes preloaded with a cryo-pack and pairs seamlessly with our pneumatic pump to deliver targeted relief. Designed to reduce pain and swelling across a wide range of procedures, this system combines smart design with proven results—making it a must-have for post-op and injury care.



Cold Compression Pro Pump comes with batteries, a screw-driver, and power supply. Cold Compression Wrap attachments are sold separately and come with one removable cryo-gel pack.

Features:

- 2 pneumatic ports to treat 1 or 2 areas at one time
- Adjustable pressure settings: 20, 40, 60, 80, or 100 mmHg to control edema
- LCD screen displays pressure to show inflation cycle
- Set treatment times; 10-minute intervals up to 120 minutes with default of 30 minutes
- Offers portable treatments with battery-powered operation, or long-duration operation via power supply—giving clinicians flexible options to support patient recovery in both mobile and stationary care settings.
- Easy-to-clean nylon material
- Removable cryo-gel pack
- Simple strapping system for donning & doffing

Suggested HCPCS: E1399

Product Name	Part Number
Cold Compression Wrap Pro Pump	CCW-PP
Cold Compression Wrap Pro Shoulder	CCW-PSW
Cold Compression Wrap Pro Knee	CCW-PKW
Cold Compression Wrap Pro Hip	CCW-PHW
Cold Compression Wrap Pro Elbow	CCW-PEW
Cold Compression Wrap Pro Back	CCW-PBW
Cold Compression Wrap Pro Ankle	CCW-PAW

For more information on suggested HCPCs and our PDAC letters, go to manamed.com/pdac-letters

2612 Sirius Drive, Denton, TX 76208

888-508-0712 (toll free)

support@Manamed.com www.Manamed.com

Manamed and the Manamed logo are trademarks of Manamed LLC.

 **MANAMED**

Recovery activated